

Getinge Clean

Universal Detergent

Mild alkaline detergent for use in washer-disinfectors, ultrasonic baths and manual cleaning of medical devices. Suitable for use on alkali resistant surfaces such as stainless steel, ceramics, glass and resistant plastics and soft metals.



Instructions for use

- Dose range 2–10 ml/lit.
- Cleaning temperature (machine) 35–60 °C.
- Cleaning temperature (manual/ultrasonic) 30–50 °C.
- Ensure complete rinsing after use.
- Neutralization is not required.

Cautions

- Do not use on polycarbonate material.
- Check material for compatibility before use.
- Do not mix with other chemicals.
- Always wear correct PPE.

Ingredients

Enzymes, sequestering agents, corrosion inhibitors and surfactants.

Storage information

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 5–25 °C. Shelf life: 18 months in original sealed containers.

Symbols



Performance tested



Detergent



All goods



Robotic surgery instruments

Technical data

Specific gravity (20 °C):	1.09 (typical)
pH as supplied (typical):	10.09
pH at 5 ml/lit water (typical):	>9.9
Viscosity (concentrate, 20 °C):	4.6 cSt

Safety and handling

Refer to Safety Data Sheet.

CE Mark

Getinge Clean Universal Detergent conforms to the European directive 93/42/EEC, Annex 1 concerning medical devices.

Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Environmental awareness

Biodegradable formulation and recyclable packaging.

Ordering number

Size	Ordering number
5 L	6001692431
205 L	6001692731

Material compatibility

Getinge Clean manual and automated detergent range		Aluminum	Anodized Aluminum	Titanium	Non-ferrous metals	Matt chrome plated	Stainless steel	Glass	Ceramic	Rubber	Plastic
■ Heavy Soil Detergent				M		•	•	•	•	•	M
■ Universal Detergent		M	M	M	M	•	•	•	•	•	M
■ Enzymatic Detergent		•	•	•	•	•	•	•	•	•	M
■ MIS Detergent		•	•	•	•	•	•	•	•	•	•
■ Prion Detergent		•	•	•	•	•	•	•	•	•	•
■ Manual Plus		•	•	•	•	•	•	•	•	•	•
■ Neutralizer		M	M	•	•	•	•	•	•	M	M
■ Neutralizer Plus				M	M	M	•	•	M	M	M
■ Chamber Foam							•				
■ Rinse Aid		•	•	•	•	•	•	•	•	M	M
■ Instrument Lubricant		•	•	•	•	•	•				
■ Instrument Lubricant Plus		•	•	•	•	•	•				
■ Pre Treatment Foam		•	•	•	•	•	•	•	•	•	•

004980

This chart is only a guideline and information on the reprocessing of medical devices from the manufacturer according to EN ISO 17664 must always be followed.

M = Always refer to instrument and washer-disinfector manufacturer's guidance.