

Getinge Clean

TEC-QUAT 128

TEC-QUAT 128 is a one-step disinfectant cleaner, sanitizer, fungicide, mildewstat, virucide* and deodorizer.

Application

This product has been formulated to aid in the reduction of cross-contamination on hard nonporous surfaces not only in hospitals, but also in institutions and industry.

Key Features

- Delivers non-acidic disinfectant and cleaning performance in an economical concentrate
- Heavy-duty disinfectant cleaner that cleans, disinfects and deodorizes in one labor-saving step
- Effective in hard water up to 200 ppm hardness in the presence of 5% serum contamination
- Formulated for use in daily maintenance programs with a balance of detergents and malodor counteractants that deliver effective cleaning, disinfecting and malodor control
- Non-scratch cleaning of showers and tubs, shower doors and curtains, fixtures and toilet bowls



Quality Statement: Confidence in Getinge is the most important quality criteria. This is the hallmark of all our external and internal commitments, activities and products. Products and services supplied by Getinge conform to the agreed terms and expectations. The achievement of these quality goals is the basis for continued competitive and successful enterprise.

Usage Information

Washable hard, nonporous surfaces:

Floors	Walls	Countertops
Cabinets	Tables	Appliances (exterior)
Desks	Doorknobs	Chairs
Tubs	Glazed Tiles	Shower Stalls
Sinks	Toilets	Bathroom Fixtures

Hard, nonporous surfaces:

Metal	Stainless Steel	Plastic Upholstery
Glazed Ceramic	Sealed Stone	Finished Woodwork
Vinyl Upholstery	Glazed Porcelain	Glazed Enamel Surfaces
Plastic (such as polystyrene or polypropylene)		

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that [1] is introduced directly to the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or [2] contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

Effective Against

This product is a one-step Hospital use disinfectant, Bactericidal according to the current AOAC Use-Dilution Test Method and Virucidal* according to the virucidal qualification modified in the presence of 200 ppm hard water plus 5% organic serum against:

*2013 Influenza A Virus	*Human Coronavirus	Staphylococcus aureus—Community Associated Methicillin-Resistant [CA-MRSA][NRS384][USA300]
*Acinetobacter baumannii	*Influenza Virus Type A / Brazil	Staphylococcus aureus—Multi-Drug Resistant [resistant to tetracycline [Tc], penicillin [Pc], streptomycin [Sm] and erythromycin [Em] and susceptible to chloramphenicol [Cm][<i>in vitro</i>]
*Adenovirus Type 5 at 1:32 dilution	*Influenza Virus Type A / Hong Kong	Staphylococcus aureus—Methicillin-Resistant (MRSA)
*Adenovirus Type 7 at 1:32 dilution	Klebsiella pneumoniae	Staphylococcus aureus—Vancomycin Intermediate Resistant (VISA)
Burkholderia cepacia	Klebsiella pneumoniae-NDM-1-positive	Streptococcus pyogenes
Campylobacter jejuni	Legionella pneumophila	Vibrio cholerae
Corynebacterium ammoniagenes	Listeria monocytogenes	*Vaccinia
Enterobacter aerogenes	*Norwalk Virus	Fungi:
Enterobacter cloacae NDM-1 positive	Pseudomonas aeruginosa	Aspergillus niger
Enterococcus faecalis [Enterococcus]	*Respiratory Syncytial Virus [RSV]	Trichophyton mentagrophytes [the Athlete's Foot Fungus][A cause of ringworm]
Enterococcus faecalis-Vancomycin resistant [VRE]	*Rotavirus	
ESBL Escherichia coli	Salmonella enterica	
Escherichia coli [E.coli]	Salmonella schottmuelleri	
Escherichia coli, New Delhi Metallo-Beta Lactamase (NDM-1)	Salmonella typhi	
*Hepatitis B Virus [HBV]	*SARS Associated Coronavirus	
*Hepatitis C Virus [HCV]	Serratia marcescens	
*Herpes Simplex Virus Type 1	Shigella dysenteriae	
*Herpes Simplex Virus Type 2	Staphylococcus aureus	
*HIV-1 [AIDS Virus]	Staphylococcus aureus—Community Associated Methicillin-Resistant [CA-MRSA][NRS123][USA400]	

***Kills HIV and HBV and HCV on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 [HIV-1] [associated with AIDS] or Hepatitis B Virus [HBV] or Hepatitis C Virus [HCV].

Instructions

Preparation of Disinfectant Use-Solution: Add 1-4 oz. per gallon of water. For heavy duty use, add 4 oz. per gallon of water [1:32]. Use 4 oz. per gallon of water to kill Adenovirus Type 5 and Type 7.

Disinfection/Virucidal*/Fungicidal/Mold and Mildew

Control Directions: Apply use-solution to hard, nonporous surfaces, thoroughly wetting surfaces with a cloth, mop, sponge, sprayer or by immersion. Treated surfaces must remain wet for 10 minutes. For Influenza Virus Type A and Human Coronavirus, treated surfaces must remain wet for 1 minute. Wipe dry with a cloth, sponge or mop, or allow to air dry. For heavily soiled areas, a preliminary cleaning is required. For sprayer applications, use a coarse spray device. Spray 6-8 inches from the surface; rub with a brush, sponge or cloth. Do not breathe spray. Rinse all surfaces that come in contact with food such as countertops, appliances, tables and stovetops with potable water before reuse. Do not use on utensils, glassware or dishes.

For fungicidal activity: This product is an effective fungicide against Trichophyton mentagrophytes when used on surfaces in areas such as locker rooms, dressing rooms, shower and bath areas and exercise facilities. Dilute at rate of 1 oz. per gallon in 200 ppm hard water and allow surface to remain wet for 10 minutes. This product is a one-step fungicide at 1 oz. per gallon in 200 ppm hard water and 5% serum.

For mold and mildew: At 1 oz. per gallon, this product will effectively inhibit the growth of mold and mildew and the odors caused by them when applied to hard, nonporous surfaces. Follow disinfection instructions. Repeat treatment every seven days, or more often if new growth appears.

To clean and deodorize: Mix 1 oz. of this product per gallon of water to clean and deodorize surfaces. Apply using a cloth, mop, sponge or sprayer. Wipe or allow to air dry.

Special Instructions

Cleaning and decontamination against HIV-1 or HBV or HCV on surfaces/objects soiled with blood/body fluids:

Personal Protection: Clean up must always be done wearing protective gloves, gowns, masks and eye protection. Cleaning Procedure: Blood and other body fluids containing HIV-1 or HBV or HCV must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Leave surface wet for 1 minute for HIV-1 and 10 minutes for HBV and HCV with a 1-4 oz. per gallon use-solution. Use 4 oz per gallon of water to kill Adenovirus type 5 and type 7. Use a 10-minute contact time for disinfection against all other viruses, bacteria and fungi claimed.

Disposal of Infectious Material: Blood, body fluids, cleaning materials and clothing must be autoclaved and disposed of according to local regulations for infectious waste disposal.

- **Do not** add bleach, detergents, or any other chemical to this product. Added chemical may destroy the effectiveness of this cleaner.

Note: Some plastics and rubbers are not compatible with detergents, disinfectants or high water temperatures. Check with device manufacturer's IFU for suitability.

Ordering Information

Getinge Clean

Order Number

4 x 1 Gallon

61301605283

5 Gallons

61301605064

Notice:

For hospital, institutional and industrial use only. Not for home use. Refer to SDS for more information regarding safety, handling, storage and waste disposal methods.



Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. • Copyright 2020 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Sales Office, US • 1 Geoffrey Way • Wayne, NJ 07470

Sales Office, Canada • 90 Matheson Blvd West, Suite 300 • Mississauga, Ontario L5R 3R3 • Canada

www.getinge.com

MCV00044188 REVC
EXP 07/22