

Suspected IAB Perforation Management

GETINGE 

1 Standby

Put pump in Standby

2 Disconnect

Disconnect helium extender tubing
from back of pump

3 Clamp

Clamp Extracorporeal Tubing
between y-fitting and male port

Clinical Considerations:

- Consider IAB catheter replacement if patient's condition warrants.
- Notify physician and prepare for possible reinsertion.
- If blood is suspected of having entered pump, take pump out of service. It should be evaluated before use in another patient by Biomed/Technical Service to determine if replacement of contaminated components are necessary.