



Intensive care ventilators

A primary option for treatment of
COVID-19 patients in the ICU

Summary and recommendations regarding ventilators

Intensive care ventilators are normally found in the hospital Intensive Care Unit. This category of ventilators supports lung-protective ventilation recommended by several international guidelines and the World Health Organization (WHO). An intensive care ventilator must be considered as the primary option during the monitoring and treatment of COVID-19 patients in the ICU.

Ventilation management functions for COVID-19	Intensive care ventilators	Emergency and transport ventilators	Solutions included in Servo-u ventilator system
Complete set of invasive ventilation modes	✓	✗	PC, PRVC, VC, PS, VS and Bi-Vent/APRV with automatic patient circuit compensation.
Complete set of non-invasive ventilation modes	✓	✗	NIV PC, NIV PS and NIV NAVA with automatic leakage compensation.
High-flow oxygen therapy	✓	✗	HFT with optional patient respiratory drive monitoring.
Basic monitoring	✓	✗	Pressure, flow and volume waveforms. Tidal volume/PBW, driving pressure, Cdyn, etCO ₂ and 72-hour, multi-parameter trending.
Advanced monitoring	✓	✗	Monitoring of patient respiratory drive (Edi, P0.1) and carbon dioxide production (VtCO ₂).
Visual decision support	✓	✗	Servo Compass including VT/PBW, plateau and driving pressure for assessment of lung-protective ventilation at a glance.
Weaning and ventilator liberation tools	✓	✗	Automode and NAVA to assure both lung- and diaphragm-protective ventilation.

Background

Mechanical ventilation is life saving for patients with respiratory failure and is a key component in the fight against the new coronavirus, severe acute respiratory syndrome coronavirus 2 (SARS-COV-2). It can also damage the lungs if the pressures or volumes are not carefully monitored and controlled.¹ Therefore, mechanical ventilation must be carefully applied in a balance between insufficient ventilation and ventilation-induced lung injury (VILI) — which has a strong impact on long-term outcome and disabilities.

Mechanical ventilation can either be invasive, via a tube in the airways, or non-invasive through a face mask or nasal prongs. However, in many cases the patient's condition will require intubation, e.g., with severe lung injury or disease.

Transport and emergency ventilators with COVID-19

When the recently mass-produced ventilators from new entrants to the industry were reviewed, it was found that most of them had a functionality comparable to today's transport or emergency ventilators (*see Appendix for classification of ventilators*). Transport and emergency ventilators are mainly used for the short term and for patients who are in stable condition. However, these ventilators unfortunately do not fulfill the technical specifications published by WHO — which describe the minimum requirements that invasive and non-invasive ventilators must comply with to ensure quality, safety and effectiveness when used for the management of COVID-19.²

It is also well-known that these simpler ventilators cannot generate the precise air delivery needed to treat patients with acute respiratory distress syndrome (ARDS)^{3,4} and they are lacking decision support tools for lung protective ventilation. In a crisis situation, these ventilators are a better alternative compared to manual ventilation ("bagging"), however they do not take the place of the intensive care ventilators.

ECRI, a worldwide organization that conducts independent medical device evaluations, clearly states that, in the case of shortages of intensive care ventilators during the COVID-19 pandemic, healthcare providers should select from available ventilation devices in the following suggested order of preference⁵:

1. Intensive care ventilators.
2. Advanced transport, sub-acute, and home care ventilators that have intensive care features and are capable of treating patients with ARDS.
3. Anesthesia machines. Anesthesia units can be re-purposed to provide ventilatory support for critically ill patients, as long as precautions are taken (*see ECRI Alert S0397: [COVID-19]*).
4. Basic transport, emergency, and home care ventilators.

Non-invasive ventilators with COVID-19

The most common diagnosis⁶ in severe COVID-19 patients is severe pneumonia associated with sepsis and ARDS⁷. The latest studies about mechanically-ventilated patients with COVID-19 show that non-invasive ventilators alone will not be sufficient for COVID-19 patients, since these type of ventilation devices are not designed for invasive ventilation of a critically ill, complex patient and most of the COVID-19 patients require invasive ventilation during their ICU stay⁸.

Intensive care ventilators with COVID-19

The first step to optimizing mechanical ventilation for patients' conditions and needs in the ICU is to use high performance ventilators⁹ such as Getinge Servo-u. The Getinge Servo-u uses built-in sensors, microprocessors, and intelligent software to adequately provide and monitor the target ventilation and automatically adjust the ventilator parameters based upon patients' needs. It is well-known that mechanical ventilation is life saving for patients with respiratory failure, but it can also damage the lungs if the pressures or volumes are not carefully controlled and monitored.

Well-recognized, key opinion leader in mechanical ventilation, Luciano Gattinoni, also focused on the importance of monitoring lung mechanics to identify the different

phenotypes of COVID-19¹⁰, and individualize the treatment accordingly. In Figure 1, drivers and interrupters of progressive lung injury in COVID-19 infection are identified.




Figure 1 Adapted from Marini JJ, Gattinoni L. Management of COVID-19 respiratory distress. JAMA. doi:10.1001/jama.2020.6825.

It is vital to know that the clinical presentation of COVID-19 patients can change over time. In their interim guidance document¹¹, WHO also strongly recommends implementing lung protective ventilation strategies for the management of critical COVID-19. Therefore, intensive care ventilators play a critical role during the COVID-19 pandemic since they:

- Can support clinicians in every phase of COVID-19, from non-invasive (NIV) to invasive ventilation modes, if the patient's condition worsens. Clinicians do not have to triage the ventilation devices at their disposal, matching the device capabilities with the severity of the patient illness due to ventilator shortages.
- Can generate the precise gas delivery needed to treat patients with Acute Lung Injury (ARDS).
- Are able to adequately provide and monitor the target ventilation and automatically adjust the ventilatory parameters, based on the patient's needs, in order to protect the lungs and prevent lung injury.
- Offer advanced monitoring tools such as P0.1, Edi. These allow the clinician to understand the effects of mechanical ventilation in real time, and adjust the ventilatory parameters on the basis of the patient's condition and needs.
- Provide early alerts and alarms to warn clinicians of a potentially harmful situation, or if the patient becomes distressed or even crashes. Please note that the alarm capabilities for emergency, transport, and home care ventilators can be very limited⁴.
- Have the ability to deliver collaborative therapy — offering nebulized drugs through a closed circuit, state-of-the-art nebulizer system.
- Allow the removal of the user interface, and keep it outside the room, in order to reduce exposure and the use of PPE.
- Have advanced weaning modalities which will help the patient liberate from mechanical ventilation faster.

Getinge’s Servo ventilators support clinicians with advanced lung protective tools — which helps to identify injurious ventilation and to diagnose clinical problems.

In Figure 2 the clinical guidelines for COVID-19, which have been developed by representatives of the University of Toronto Interdepartmental Division of Critical Care Medicine, are demonstrated¹². Using these guidelines we are able to show where Getinge can support or offer a solution.

COVID-19 Initial Management

How severe is the hypoxemia?
(PaO_2/FiO_2 or $SpO_2/FiO_2 >150$ mmHg)

Oxygen therapy via face mask or High Flow Therapy.
Target SpO_2 : 90–94%



Worsening oxygenation and other indications for intubation →

Treatment Algorithm

Intubation and controlled mechanical ventilation

Ventilator settings and targets:
Tidal Volume: 6–8 ml/kg PBW
Plateau pressure <30 cmH₂O
Driving pressure <15cmH₂O
PEEP: 8–10 cmH₂O
Maintain deep sedation and consider NMB



Assess for recruitability

No →

Consider prone positioning



↓ Yes



Recruitment maneuver

↓ Oxygenation not getting better
($PaO_2/FiO_2 \leq 150$ mmHg)



Consider prone positioning

↓ If refractory hypoxemia

→

Consider ECLS



Getinge has the functionality to support

Figure 2

References

1. Dreyfuss D, Basset G, Soler P, Saumon G. Intermittent positive-pressure hyperventilation with high inflation pressures produces pulmonary microvascular injury in rats. *Am Rev Respir Dis*. 1985 Oct;132(4):880-4.
2. Technical specifications for invasive and non-invasive ventilators for COVID-19: interim guidance: https://apps.who.int/iris/bitstream/handle/10665/331792/WHO-2019-nCoV-Clinical-Ventilator_Specs-2020-1-eng.pdf?sequence=1&isAllowed=y.
3. Chipman DW, Caramenz MP, Miyoshi E, et al. Performance comparison of 15 transport ventilators. *Respir Care*. 2007 Jun;52(6):740-51.
4. ECRI Exclusive Special Report - S0404: [COVID-19] Philips— Respiroics E30 Ventilators: ECRI Assessment of Emergency Use Authorization Device.
5. ECRI Exclusive Special Report - S0398 : [COVID-19] Shortages of Intensive Care Ventilators—Strategies for Mitigation. Medical Device Special Report.
6. Clinical management of severe acute respiratory infection (SARI) when COVID-19 disease is suspected Interim guidance 13 March 2020. WHO reference number: WHO/2019-nCoV/clinical/2020.4.
7. Zhu J, Ji P, Pang J, et al. Clinical characteristics of 3062 COVID-19 patients: A meta-analysis. *J Med Virol*. 2020 Apr 15;10.1002/jmv.25884.
8. Ziehr DR, Alladina J, Petri CR, et al. Respiratory Pathophysiology of Mechanically Ventilated Patients with COVID-19: A Cohort Study. *Am J Respir Crit Care Med*. 2020 Jun 15;201(12):1560-1564.
9. Govoni L, Dellaca RL, Peñuelas O, et al. Actual performance of mechanical ventilators in ICU: a multicentric quality control study. *Med Devices (Auckl)*. 2012;5:111-9.
10. Gattinoni L, Chiumello D, Caironi P, et al. COVID-19 pneumonia: different respiratory treatment for different phenotypes? *Intensive Care Med*. 2020 Jun;46(6):1099-1102.
11. Analysis of the technologies required for clinical management of COVID-19 patients <https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/patient-management>.
12. <https://www.criticalcare.utoronto.ca/covid-19-resources>.

Appendix: Classification of ventilators

- **Intensive care ventilators:** Ventilators designed for use in the ICU.
- **Home care ventilators:** Ventilators used for long-term ventilation in the patient's home environment.
- **Sub-acute ventilators:** Ventilators that are designed for use outside of the ICU, through either non-invasive ventilation, or ventilation of stable patients. This includes step-down wards, general wards, and skilled nursing facilities.
- **Emergency & Transport ventilators:** Ventilators used in emergency situations both within the hospital and externally (mainly for short-term).



Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Automode, NAVA, Servo-u and Servo Compass are trademarks of Maquet Critical Care AB • Copyright 2020 Getinge AB or its subsidiaries or affiliates • All rights reserved • ⚠ CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician. Refer to Instructions for Use for current indications, warnings, contraindications, and precautions.

Sales Office • Getinge • 1 Geoffrey Way • Wayne, NJ 07470 • USA

Manufacturer • Maquet Critical Care AB • 171 54 Solna, Sweden • Phone: +46 (0)10 335 00 00 • info@getinge.com

www.getinge.com

MX-7818
MCV00105929 REV B
EXP08/22