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
**Lucea 50-100**  
**Lucea 50**

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| 국가                  | 참조번호                             | 버전   | 제목                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|---------------------|----------------------------------|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| EU                  | ROHS 지침                          | 2011 | DIRECTIVE 2011/65/EU OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment                                                                                                                                                                                                                                                         |
|                     |                                  | 2015 | COMMISSION DELEGATED DIRECTIVE (EU) 2015/863 of 31 March 2015, amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances                                                                                                                                                                                                                                            |
|                     |                                  | 2016 | COMMISSION DELEGATED DIRECTIVE (EU) 2016/585 of 12 February 2016 amending, for the purposes of adapting to technical progress, Annex IV to Directive 2011/65/EU of the European Parliament and of the Council as regards an exemption for lead, cadmium, hexavalent chromium, and polybrominated diphenyl ethers (PBDE) in spare parts recovered from and used for the repair or refurbishment of medical devices or electron microscopes      |
|                     |                                  | 2017 | DIRECTIVE (EU) 2017/2102 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 15 November 2017 amending Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment                                                                                                                                                                                                                  |
| 전 세계                | IEC 63000                        | 2022 | Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances                                                                                                                                                                                                                                                                                                       |
| EU                  | REACH Regulation                 | 2006 | REGULATION (EC) No 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 18 December 2006 concerning the Registration, Evaluation, Authorisation and REACH - Restriction of Chemicals (REACH), amending Directive 1999/45/EC and repealing Council Regulation (EEC) No 793/93 and Commission Regulation (EC) No 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC |
| USA _<br>California | US California proposition 65 Act | 1986 | HEALTH AND SAFETY CODE - HSC DIVISION 20. MISCELLANEOUS HEALTH AND SAFETY PROVISIONS CHAPTER 6.6. Safe Drinking Water and Toxic Enforcement Act of 1986                                                                                                                                                                                                                                                                                        |
| 중국                  | SJ/T 11365-2006                  | 2006 | ACPEIP - Administrative Measure on the Control of Pollution caused by Electronic Information Products Chines RoHS (Restriction of Hazardous Substances)                                                                                                                                                                                                                                                                                        |

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| 국가                     | 참조번호                      | 연도   | 제목                                                                                                                                                                                                                                                               |
|------------------------|---------------------------|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 아르헨티나                  | Dispocision 2318/2002     | 2002 | Administración Nacional de Medicamentos, Alimentos y Tecnología Médica - Registro de productos Medicas - Reglamento                                                                                                                                              |
| 호주                     | TGA 236-2002              | 2021 | Therapeutic Goods (Medical Devices) Regulations 2002. Statutory Rules No. 236, 2002 made under the Therapeutic Goods Act 1989                                                                                                                                    |
| Bosnia and Herzegovina | Act                       | 2008 | Medicinal products and medical devices act of Bosnia and Herzegovina ("Official Gazette of BiH, No. 58/08)                                                                                                                                                       |
| 브라질                    | RDC 665/2022              | 2022 | Resolution RDC n°665, 30 March 2022, Provides for the Good Manufacturing Practices for Medical Devices and Medical devices for In Vitro Diagnostis                                                                                                               |
| 브라질                    | RDC 751/2022              | 2022 | RDC No. 751, of September 15, 2022, which provides for risk classification, notification and registration regimes, and labeling requirements and instructions for use of medical devices.                                                                        |
| 브라질                    | Ordinance 384/2020        | 2020 | INMETRO Certification - Compliance Assessment Requirements for Equipment under Health Surveillance Regimen - Consolidated.                                                                                                                                       |
| 캐나다                    | SOR/98-282                | 2024 | Medical Devices Regulations                                                                                                                                                                                                                                      |
| 중국                     | Regulation n°739          | 2021 | Regulation for the Supervision and Administration of Medical Devices                                                                                                                                                                                             |
| 콜롬비아                   | Decree 4725               | 2005 | DECRETO NÚMERO 4725 DE 2005 (Diciembre 26) por el cual se reglamenta el régimen de registros sanitarios, permiso de comercialización y vigilancia sanitaria de los dispositivos médicos para uso humano.                                                         |
| EU                     | Regulation 2017/745/EU    | 2017 | REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC |
| 인도                     | Rule                      | 2017 | Medical Device Rules, 2017                                                                                                                                                                                                                                       |
| 인도네시아                  | Regulation 62             | 2017 | Regulation of the minister of health of the republic of Indonesia number 62 of 2017 on product license of medical devices, in vitro diagnostic medical devices and household health products                                                                     |
| 이스라엘                   | Law 5772-2012             | 2012 | The Medical Equipment Law, 5772-2012                                                                                                                                                                                                                             |
| 일본                     | MHLW Ordinance: MO n °169 | 2021 | Ministerial Ordinance on Standards for Manufacturing Control and Quality Control for Medical Devices and In-Vitro Diagnostics                                                                                                                                    |
| 케냐                     | Act                       | 2002 | The Pharmacy and Poisons Act, Cap 244 of the Laws of Kenya                                                                                                                                                                                                       |
| 말레이시아                  | Act 737                   | 2012 | Medical Device Act 2012 (Act 737)                                                                                                                                                                                                                                |
| Montenegro             | Law 53/09                 | 2009 | Law of Montenegro on Medical Devices (2009)                                                                                                                                                                                                                      |
| 모로코                    | Law 84-12                 | 2012 | Law n°84-12 relative to medical devices                                                                                                                                                                                                                          |

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| 국가       | 참조번호                | 연도   | 제목                                                                                                                                                                                                                                        |
|----------|---------------------|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 뉴질랜드     | Regulation 2003/325 | 2003 | Medicines (Database of Medical Devices) Regulations 2003 (SR 2003/325)                                                                                                                                                                    |
| 사우디 아라비아 | Regulation          | 2017 | "Medical Device Interim Regulation" issued by the Board of Directors of the Food and Drug Authority (1-8-1429) dated 29/12/1429 H and amended by Saudi Food and Drug Authority Board of Directors decree No. (4-16-1439) dated 27/12/2017 |
| 세르비아     | Law 105/2017        | 2017 | Law on Medicinal Products and Medical Devices, "Official Gazette of the Republic of Serbia," No. 105/2017                                                                                                                                 |
| 대한민국     | 법률 제14330호          | 2016 | 의료기기법                                                                                                                                                                                                                                     |
| 대한민국     | 법령 제27209호          | 2016 | 의료법 시행령                                                                                                                                                                                                                                   |
| 대한민국     | 규칙 제1354호           | 2017 | 의료법 시행규칙                                                                                                                                                                                                                                  |
| 스위스      | RS (Odim) 812.213   | 2020 | Medical Devices Ordinance (MedDO) of 1 July 2020                                                                                                                                                                                          |
| 대만       | TPAA 2018-01-31     | 2018 | Taiwanese Pharmaceutical Affairs Act                                                                                                                                                                                                      |
| 태국       | Act 2562            | 2019 | Medical Device Act (No. 2) B.E. 2562(2019)                                                                                                                                                                                                |
| 영국       | Act                 | 2021 | Medical Devices Regulations 2002 n°618                                                                                                                                                                                                    |
| 미국       | 21CFR Part 7        | 2023 | Title 21--Food And Drugs<br>Chapter I--Food And Drug Administration<br>Department Of Health And Human Services<br>Subchapter A -- General<br>PART 7 - Enforcement policy                                                                  |
| 미국       | 21CFR Subchapter H  | -    | Title 21--Food And Drugs<br>Chapter I--Food And Drug Administration<br>Department Of Health And Human Services<br>Subchapter H -- Medical Devices                                                                                         |
| 베트남      | Decree 98/2021      | 2021 | Decree No. 98/2021/ND-CP November 8, 2021 of the Government on the management of medical equipment                                                                                                                                        |

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© Copyright 2024
Maquet SAS · Parc de Limère · Avenue de la Pomme de Pin · CS 10008 ARDON ·