Cardiosave Hybrid / Rescue IABP

Instructions For Use Addendum

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1. Cardiosave Hybrid & Rescue Operating Instructions Addendum

1.1 Schedule B Update:

The following Schedule B table entries have been added and updated:

Schedule B

Required Action (Refer to Service Manual)	Every 12 Months or 2500 Hours*	6,000,000 Inflate, Deflate Cycles or 4 Years.*	12,000,000 Inflate, Deflate Cycles, Every 2000 Assist Hours, or 4 Years.*	Every 5000 Hours	Every 8 Years
Replace the Real Time Clock NVRAM IC on the Executive Processor Board.					●
Replace the quick disconnect fitting o-ring and lubricate the console release latch.	●				

* Whichever comes first.

1.2 Battery Pack Run Time Specification Update:

The following Battery Pack Run Time specification has been updated and note message has been added:

Battery Pack (Cardiosave supports the use of two battery packs)

(P/N 0146-00-0097)

Run Time:	60 Minutes Minimum, 120 bpm, 22 +/- 5 degrees C
(per pack, with full charge)	80 Minutes Typical - new battery, 120 bpm, 22 +/- 5 degrees C

NOTE

For a new battery, 80 Minutes at 120 bpm is equivalent to 90 Minutes at 90 bpm per previous battery run time claims.

NOTE

A reduction in run time can occur over a battery's life for reasons such as age, storage temperature and discharge depth. Environmental conditions and sustained heart rates above 120 bpm may affect run times. Monitor battery charge levels closely during these conditions.

The 0146-00-0097 **Cardiosave** Lithium-Ion Battery Pack is distinguished by the part number revision and corresponding serial number suffix.

Revision P or earlier: Identified by a serial number suffix of "**IP**": Nominal voltage of 15.0 VDC, Capacity of 14.3 Ah.

Revision Q or later: Identified by a serial number suffix of "**ACCUTR**": Nominal voltage of 14.4 VDC, Capacity of 16.8 Ah.

Battery Pack (0146-00-0097) Weight: 1.4 kg $\pm$ 5%, 3.1 lbs $\pm$ 5%

NOTE

Both versions are forwards and backwards compatible with the **Cardiosave** IABP system.

1.3

Vibration and Shock Table Specifications Update:

The following Vibration and Shock specifications have been updated:

1.3.1

Vibration/Shock (2nd Edition)

Sinusoidal Vibration:	EN 60068-2-6:2008 - Environmental Testing, Part 2 Test Fc: Vibration (sinusoidal) EN13718-1:2008 (E) section 4.6.2 which cites ISO 7137:1995. Test method is defined per RTCA/DO-160F par 8.5.1 Cat S, fuselage zone M.
Random Vibration:	IEC 60068-2-64:2019, frequency range 20 to 500 Hz, ASD $0.02g^2/Hz$, duration 9 min. Nonoperational. For ground applications: EN 1789:2007 "Medical vehicles and their equipment - Road ambulances". For aircraft transport: RTCA/DO-160F "Environmental Conditions and Test Procedures for Airborne Equipment".
Bump Test:	EN 1789:2007 - Medical Vehicles and their equipment - Road ambulances: Functional shock test per IEC 60068-2-29, Test Eb, 15g, duration 6 ms, 1000 times in normal operating position.
Rough Handling Test:	AAMI: 60601-1:2005
Shock - Operational:	RTCA/DO-160F Shock, Category B Operational standard test for fixed wing and helicopter: 6g, duration 11ms, sawtooth, 3 times in each orientation
Shock - Non-Operational:	EN60068-2-27:2008, Environmental Testing, Part 2 Test Ea and guidance: Shock per Table 1: 20g, duration 11 ms, saw-tooth, one in each orientation (6 total, can use dummy load) EN 1789:2007 - Medical Vehicles and their Equipment section 6.3.5.
Crash Test:	RTCA/DO-160F - Environmental Conditions and Test Procedures for Airborne Equipment: Crash safety test for fixed wing and helicopter, 20g, duration 6ms, saw-tooth, one in each orientation (6 total, can use dummy load)
Elevator Threshold Test:	AAMI: 60601-1:2005
Barrier Impact Test:	AAMI: 60601-1:2005
Drop Test:	AAMI: 60601-1:2005
Tip-over Test:	UL 60601-1:2006

1.3.2

Vibration/Shock (3rd Edition)

Sinusoidal Vibration:	EN 60068-2-6:2008 - Environmental Testing, Part 2 Test Fc: Vibration (sinusoidal) EN13718-1:2008 (E) section 4.6.2 which cites ISO 7137:1995. Test method is defined per RTCA/DO-160F par 8.5.1 Cat S, fuselage zone M.
Random Vibration:	IEC 60068-2-64:2019, frequency range 20 to 500 Hz, ASD $0.02g^2/Hz$, duration 9 min. Nonoperational. For ground applications: EN 1789:2007 "Medical vehicles and their equipment - Road ambulances". For aircraft transport: RTCA/DO-160F "Environmental Conditions and Test Procedures for Airborne Equipment".
Bump Test:	EN 1789:2007 - Medical Vehicles and their equipment - Road ambulances: Functional shock test per IEC 60068-2-27:2008, Test Eb, 15g, duration 6 ms, 1000 times in normal operating position.
Rough Handling Test:	IEC 60601-1:2012
Shock - Operational:	RTCA/DO-160F Shock, Category B Operational standard test for fixed wing and helicopter: 6g, duration 11ms, sawtooth, 3 times in each orientation
Shock - Non-Operational:	EN60068-2-27:2008, Environmental Testing, Part 2 Test Ea and guidance: Shock per Table 1: 20g, duration 11 ms, saw-tooth, one in each orientation (6 total, can use dummy load) EN 1789:2007 - Medical Vehicles and their Equipment section 6.3.5.
Crash Test:	RTCA/DO-160F - Environmental Conditions and Test Procedures for Airborne Equipment: Crash safety test for fixed wing and helicopter, 20g, duration 6ms, saw-tooth, one in each orientation (6 total, can use dummy load)
Elevator Threshold Test:	IEC 60601-1:2012
Barrier Impact Test:	IEC 60601-1:2012
Drop Test:	IEC 60601-1:2012
Tip-over Test:	IEC 60601-1:2012

NOTE

Cardiosave IABP Hardware starting with serial number prefix CH are 60601-1 2nd Edition compliant, prefixes CB, CM and CD are 60601-1 3rd Edition compliant. The **Cardiosave** Software is 60601-1 3rd Edition compliant starting with revisions C.03 and above. Both the IABP hardware and software must be 3rd edition compliant for the **Cardiosave** IABP to be 60601-1 3rd Edition compliant.

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US

Getinge is a global provider of innovative solutions for operating rooms, intensive care units, sterilization departments and for life science companies and institutions. Based on our firsthand experience and close partnerships with clinical experts, healthcare professionals and medtech specialists, we are improving the everyday life for people - today and tomorrow.

Datascope Corp. · 1300 MacArthur Blvd., Mahwah, NJ 07430, USA · Phone: +1 800 777 4222 or +1 201 995 8700

PN: MCA000000751 Rev B · GETINGE ✱, Getinge, Maquet, **Cardiosave**, **Cardiosave Hybrid**, and **Cardiosave Rescue** are trademarks or registered trademarks of Getinge AB, its subsidiaries, or affiliates in the United States or other countries. · **Cardiosave** and Maquet are registered with the U.S. Patent and Trademark Office. · All other company or product names are the trademarks or registered trademarks of their respective holders. · Copyright 2025 Datascope Corp. · All rights reserved. · Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician. · ⚠ Refer to Instructions for Use for current indications, warnings, contraindications, and precautions. · 07/2025

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All rights reserved. No part of this publication may be duplicated, adapted or translated without prior written permission, except within the framework of the copyright laws.
© Copyright 2025 Datascope Corp.
Datascope Corp. · 1300 MacArthur Blvd., Mahwah, NJ 07430, VS · Telefoon: +1 800 777 4222 of +1 201 995 8700
Datascope Corp. · 1300 MacArthur Blvd., Mahwah, NJ 07430, USA · Telefon: +1 800 777 4222 eller +1 201 995 8700
Datascope Corp. · 1300 MacArthur Blvd., Mahwah, NJ 07430, USA · Tel: +1 800 777 4222 lub +1 201 995 8700
Datascope Corp. · 1300 MacArthur Blvd., Mahwah, NJ 07430, USA · Telefon: +1 800 777 4222 vagy +1 201 995 8700