

Getinge Clean

TEC-QUAT 128

TEC-QUAT 128 is a one-step disinfectant cleaner, sanitizer, fungicide, mildewstat, virucide* and deodorizer.



Application

This product has been formulated to aid in the reduction of cross-contamination on hard nonporous surfaces not only in hospitals, institutions and industry.

Key Features

- Delivers non-acidic disinfectant and cleaning performance in an economical concentrate
- Heavy-duty disinfectant cleaner that cleans, disinfects and deodorizes in one labor-saving step
- Effective in hard water up to 200 ppm hardness in the presence of 5% serum contamination
- Formulated for use in daily maintenance programs with a balance of detergents, biocides and malodor counteractants that deliver effective cleaning, disinfecting and malodor control
- Non-scratch cleaning of showers and tubs, shower doors and curtains, fixtures and toilet bowls

Quality Statement: Confidence in Getinge is the most important quality criteria. This is the hallmark of all our external and internal commitments, activities and products. Products and services supplied by Getinge conform to the agreed terms and expectations. The achievement of these quality goals is the basis for continued competitive and successful enterprise.

Usage Information

Washable hard, nonporous surfaces:

| | | |
|----------|--------------|-----------------------|
| Floors | Walls | Countertops |
| Cabinets | Tables | Appliances (exterior) |
| Desks | Doorknobs | Chairs |
| Tubs | Glazed Tiles | Shower Stalls |
| Sinks | Toilets | Bathroom Fixtures |

Hard, nonporous surfaces:

| | | |
|--|------------------|------------------------|
| Metal | Stainless Steel | Plastic Upholstery |
| Glazed Ceramic | Sealed Stone | Finished Woodwork |
| Vinyl Upholstery | Glazed Porcelain | Glazed Enamel Surfaces |
| Plastic (such as polystyrene or polypropylene) | | |

Effective Against

This product is a one-step Hospital use disinfectant, Bactericidal according to the current AOAC Use-Dilution Test Method and Virucidal* according to the virucidal qualification modified in the presence of 200 ppm hard water plus 5% organic serum against:

- *2013 Influenza A Virus
- *Acinetobacter baumannii
- *Adenovirus Type 5 at 1:32 dilution
- *Adenovirus Type 7 at 1:32 dilution
- Burkholderia cepacia
- Campylobacter jejuni
- Corynebacterium ammoniagenes
- Enterobacter aerogenes
- Enterobacter cloacae NDM-1 positive
- Enterococcus faecalis [Enterococcus]
- Enterococcus faecalis- Vancomycin resistant [VRE]
- ESBL Escherichia coli
- Escherichia coli [E.coli]
- Escherichia coli, New Delhi Metallo-Beta Lactamase (NDM-1)
- *Hepatitis B Virus [HBV]
- *Hepatitis C Virus [HCV]
- *Herpes Simplex Virus Type 1
- *Herpes Simplex Virus Type 2
- *HIV-1 [AIDS Virus]
- *Human Coronavirus
- *Influenza Virus Type A / Brazil
- *Influenza Virus Type A / Hong Kong
- Klebsiella pneumoniae
- Klebsiella pneumoniae-NDM-1-positive
- Legionella pneumophila
- Listeria monocytogenes
- *Norwalk Virus
- Pseudomonas aeruginosa
- *Respiratory Syncytial Virus [RSV]
- *Rotavirus
- Salmonella enterica
- Salmonella schottmuelleri
- Salmonella typhi

- *SARS Associated Coronavirus
- Serratia marcescens
- Shigella dysenteriae
- Staphylococcus aureus
- Staphylococcus aureus—Community Associated Methicillin-Resistant [CA-MRSA][NRS123][USA400]
- Staphylococcus aureus—Community Associated Methicillin-Resistant [CA-MRSA][NRS384][USA300]
- Staphylococcus aureus—Multi-Drug Resistant [resistant to tetracycline [Tc], penicillin [Pc], streptomycin [Sm] and erythromycin [Em] and susceptible to chloramphenicol [Cm] in vitro]
- Staphylococcus aureus—Methicillin-Resistant (MRSA)
- Staphylococcus aureus—Vancomycin Intermediate Resistant (VISA)
- Streptococcus pyogenes
- Vibrio cholerae
- *Vaccinia
- Fungi:
 - Aspergillus niger
 - Trichophyton mentagrophytes [the Athlete's Foot Fungus][A cause of ringworm]

Usage Information

Preparation of Disinfectant Use-Solution: Add 1-4 oz. per gallon of water. For heavy duty use, add 4 oz. per gallon of water [1:32]. Use 4 oz. per gallon of water to kill Adenovirus Type 5 and type 7.

Disinfection/Virucidal*/Fungicidal/Mold and Mildew

Control Directions: Apply use-solution to hard, nonporous surfaces, thoroughly wetting surfaces with a cloth, mop, sponge, sprayer or by immersion. Treated surfaces must remain wet for 10 minutes. For Influenza Virus Type A and Human Coronavirus, treated surfaces must remain wet for 1 minute. Wipe dry with a cloth, sponge or mop, or allow to air dry. For heavily soiled areas, a preliminary cleaning is required. For sprayer applications, use a coarse spray device. Spray 6-8 inches from the surface; rub with a brush, sponge or cloth. Do not breathe spray. Rinse all surfaces that come in contact with food such as countertops, appliances, tables and stovetops with potable water before reuse. Do not use on utensils, glassware or dishes. Prepare a fresh solution daily or more often if the solution becomes visibly dirty or diluted.

For fungicidal activity: This product is an effective fungicide against Trichophyton mentagrophytes when used on surfaces in areas such as locker rooms, dressing rooms, shower and bath areas and exercise facilities. Dilute at rate of 1 oz. per gallon in 200 ppm hard water and allow surface to remain wet for 10 minutes. This product is a one-step fungicide at 1 oz. per gallon in 20ppm hard water and 5% serum.

For mold and mildew: At 1 oz. per gallon, this product will effectively inhibit the growth of mold and mildew and the odors caused by them when applied to hard, nonporous surfaces. Follow disinfection instructions. Repeat treatment every seven days, or more often if new growth appears.

This product is not to be used as a terminal sterilant/ high level disinfectant on any surface or instrument that [1] is introduced directly to the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or [2] contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

To clean and deodorize: Mix 1 oz. of this product per gallon of water to clean and deodorize surfaces. Apply using a cloth, mop, sponge or sprayer. Wipe or allow to air dry.

***Kills HIV and HBV and HCV on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 [HIV-1] [associated with AIDS] or Hepatitis B Virus [HBV] or Hepatitis C Virus [HCV].

Special Instructions for Cleaning and Decontamination against HIV-1 or HBV or HCV on Surfaces/Objects Soiled with Blood/Body Fluids

Personal Protection: Clean up must always be done wearing protective gloves, gowns, masks and eye protection.

Cleaning Procedure: Blood and other body fluids containing HIV-1 or HBV or HCV must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Leave surface wet for 1 minute for HIV-1 and 10 minutes for HBV and HCV with a 1-4 oz. per gallon use-solution. Use 4 oz per gallon of water to kill Adenovirus type 5 and type 7. Use a 10-minute contact time for disinfection against all other viruses, bacteria and fungi claimed.

Disposal of Infectious Material: Blood, body fluids, cleaning materials and clothing must be autoclaved and disposed of according to local regulations for infectious waste disposal.

- Do not add bleach, detergents, or any other chemical to this product. Added chemical may destroy the effectiveness of this cleaner.

Note: Some plastics and rubbers are not compatible with detergents, disinfectants or high water temperatures. Check with device manufacturer's IFU for suitability.

Notice: For hospital, institutional and industrial use only. Not for home use. Refer to SDS for more information regarding safety, handling, storage and waste disposal methods.

Ordering Information

Getinge Clean

Order Number

4 x 1 Gallon

61301605283

5 Gallons

61301605064

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2017 Getinge AB or its subsidiaries or affiliates • All rights reserved. • MCV00044188 REVA

Manufacturer: Intercon Chemical Co, 1100 Central Industrial Dr., St. Louis, MO 63110
Sales Office: Getinge, 1777 East Henrietta Road, Rochester, NY 14623

www.getinge.com

GETINGE 