

Getinge Clean Chamber Foam

Foam based, ready to use acid cleaner. Designed to remove scale, mineral deposits and rust from stainless steel surfaces of autoclaves and automated washer chambers. Can also be used on stainless steel loading equipment or any other stainless steel surface. Restores stainless steel surface finish.



Instructions for use

Wear correct PPE. Apply foam to surface, leaving for up to 30 min; if necessary brush (acid resistant, non-abrasive brush). Repeat if required. Remove all foam residues by rinsing with water and/or wet cloths. Run a process without load.

Cautions

- Always wear correct PPE
- Do not apply to hot surfaces
- Rinse thoroughly
- Intended for use on stainless steel only
- Autoclaves: Avoid product entering drain

Symbols



Ingredients

Water, phosphoric acid, nitric acid, rheology modifier, surfactants and corrosion inhibitor.

Technical data

Specific gravity (20 °C):	1.024
pH as supplied (typical):	1-2

Storage information

Store in a cool, dry place protected from frost and away from alkalis and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 5–25 °C. Shelf life: 24 months in original sealed containers.

Safety and handling

Refer to Safety Data Sheet.

Environmental awareness

Biodegradable formulation and recyclable packaging.

CE Mark

Getinge Clean Chamber Foam conforms to the European directive 93/42/EEC, Annex 1 concerning medical devices.

Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Ordering number

Size	Language	Ordering number	
500 ml	EN	6002000801	To be ordered in multiples of 6

Material compatibility

Getinge Clean manual and automated detergent range		Material									
		Aluminum	Anodized Aluminum	Titanium	Non-ferrous metals	Matt chrome plated	Stainless steel	Glass	Ceramic	Rubber	Plastic
■ Heavy Soil				M		•	•	•	•	•	M
■ Universal		M	M	M	M	•	•	•	•	•	M
■ Enzymatic		•	•	•	•	•	•	•	•	•	M
■ MIS		•	•	•	•	•	•	•	•	•	•
■ Prion		•	•	•	•	•	•	•	•	•	•
■ Manual Plus		•	•	•	•	•	•	•	•	•	•
■ Neutralizer		M	M	•	•	•	•	•	•	M	M
■ Neutralizer Plus				M	M	M	•	•	M	M	M
■ Chamber Foam							•				
■ Rinse Aid		•	•	•	•	•	•	•	•	M	M
■ Instrument Lubricant		•	•	•	•	•	•				
■ Instrument Lubricant Plus		•	•	•	•	•	•				
■ Pre-treatment Foam		•	•	•	•	•	•	•	•	•	•

004980

This chart is only a guideline and information on the reprocessing of medical devices from the manufacturer according to EN ISO 17664 must always be followed.

M = Always refer to instrument and washer-disinfector manufacturer's guidance.