

## GETINGE CLEAN

### Product Specification

### Pre Treatment Foam

#### PRODUCT DESCRIPTION

A ready to use pre-treatment, pH neutral waterbased foam to be applied on medical devices immediately after use to prevent drying of surgical soils. Adheres to medical device surfaces and keeps soils moist for easy removal.

Suitable for stainless steel, ceramics, glass, plastic and soft metal surfaces. Water based and easily rinsed free in automated or manual pre-wash applications.

#### DIRECTION OF USE

Spray on soiled instruments at the immediately after use prior to transport.

Please refer to product safety data sheet before use. Use only as directed.

Do not mix with other chemicals. Shake before use.

#### INGREDIENTS

Surfactants, rheology modifier, preservative and corrosion inhibitors.

#### TECHNICAL DATA

Specific gravity (20 °C): 1.016

pH as supplied (typical): 7

pH at 5ml/water (typical): N/A

Viscosity: 110 cps (Typical)

#### STORAGE INFORMATION

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents.

Store upright in original containers.

Optimum storage temperature 5-25°C.

#### SAFETY AND HANDLING

Refer to Material Safety Data Sheet.

#### CE MARK

Instrument pretreatment conforms to the European directive 93/42/EC, Annex 1 concerning medical devices.

#### ENVIRONMENTAL AWARENESS

Biodegradable formulation and recyclable packaging.



Instruments



Foam



Getinge's range of detergents have been developed, tested and validated on Getinge washers to consistently meet ISO 15883 standards for cleaning efficacy of medical devices.



## GETINGE CLEAN

Product Specification

### Pre Treatment Foam

#### SERVICE

For technical support in regards to appropriate detergent dosing and optimization of your cleaning processes, please contact your local Getinge representative.

#### GETINGE MANUAL & AUTOMATED DETERGENT RANGE



|                                      | Aluminum | Anodised aluminium | Titanium | Non ferrous metals | Matt chrome plated | Stainless steel | Glass | Ceramic | Rubber | Plastic |
|--------------------------------------|----------|--------------------|----------|--------------------|--------------------|-----------------|-------|---------|--------|---------|
| ■ GETINGE CLEAN HEAVY SOIL DETERGENT |          |                    | M        |                    | ●                  | ●               | ●     | ●       | ●      | M       |
| ■ GETINGE CLEAN UNIVERSAL DETERGENT  | M        | M                  | M        | M                  | ●                  | ●               | ●     | ●       | ●      | M       |
| ■ GETINGE CLEAN ENZYMATIC DETERGENT  | ●        | ●                  | ●        | ●                  | ●                  | ●               | ●     | ●       | ●      | M       |
| ■ GETINGE CLEAN NEUTRALIZER          | M        | M                  | ●        | ●                  | ●                  | ●               | ●     | ●       | M      | M       |
| ■ GETINGE CLEAN RINSE AID            | ●        | ●                  | ●        | ●                  | ●                  | ●               | ●     | ●       | M      | M       |
| ■ GETINGE CLEAN INSTRUMENT LUBRICANT | ●        | ●                  | ●        | ●                  | ●                  | ●               |       |         |        |         |
| ■ GETINGE CLEAN PRE TREATMENT FOAM   | ●        | ●                  | ●        | ●                  | ●                  | ●               | ●     | ●       | ●      | ●       |

This chart is only a guideline and information on the reprocessing of medical devices from the manufacturer according to EN ISO 17664 must always be followed.

M=Always refer to instrument manufacturer's guidance.

#### ORDERING NUMBER

|           |
|-----------|
| 500 ml    |
| 503447202 |

# GETINGE

GETINGE GROUP

Getinge Disinfection AB, Ljungadalsgatan 11  
 Box 1505, SE-351 15 Växjö, Sweden  
 Phone: +46 (0)10 335 98 00  
 Web: <http://www.getinge.com>  
 E-mail: [info.disinfection@getinge.com](mailto:info.disinfection@getinge.com)