

# Ordering information


## Advanta V12 Balloon Expandable Covered Stent

5–10 mm diameter, .035" guidewire

Stent diameter / length	Order number 80 cm catheter length	Order number 120 cm catheter length	Foreshortened length		Introducer compatibility
			8 ATM nominal pressure	12 ATM rated burst pressure	
5 x 16 mm	85340	85350	15.9 mm	15.6 mm	6 Fr
5 x 22 mm	85341	85351	21.3 mm	21.0 mm	6 Fr
5 x 32 mm	85388	85394	32.3 mm	32.3 mm	7 Fr
5 x 38 mm	85320	85330	37.2 mm	37.7 mm	7 Fr
5 x 59 mm	85321	85331	58.6 mm	60.0 mm	7 Fr
6 x 16 mm	85342	85352	15.7 mm	15.1 mm	6 Fr
6 x 22 mm	85343	85353	20.8 mm	20.2 mm	6 Fr
6 x 32 mm	85389	85395	31.7 mm	31.5 mm	7 Fr
6 x 38 mm	85322	85332	36.6 mm	37.0 mm	7 Fr
6 x 59 mm	85323	85333	57.8 mm	58.7 mm	7 Fr
7 x 16 mm	85344	85354	15.0 mm	14.2 mm	7 Fr
7 x 22 mm	85345	85355	20.1 mm	19.4 mm	7 Fr
7 x 32 mm	85390	85396	31.3 mm	31.2 mm	7 Fr
7 x 38 mm	85324	85334	35.8 mm	35.7 mm	7 Fr
7 x 59 mm	85325	85335	57.1 mm	57.5 mm	7 Fr
8 x 32 mm	85391	85397	30.0 mm	29.6 mm	7 Fr
8 x 38 mm	85326	85336	34.7 mm	34.7 mm	7 Fr
8 x 59 mm	85327	85337	56.0 mm	56.5 mm	7 Fr
9 x 32 mm	85392	85398	28.7 mm	29.2 mm	7 Fr
9 x 38 mm	85328	85338	33.7 mm	32.7 mm	7 Fr
9 x 59 mm	85329	85339	54.6 mm	54.0 mm	7 Fr
10 x 38 mm	85360	85364	30.8 mm	30.9 mm	7 Fr
10 x 59 mm	85361	85365	53.3 mm	52.5 mm	7 Fr

12 mm large diameter, .035" guidewire

Stent diameter / length	Order number 80 cm catheter length	Order number 120 cm catheter length	Foreshortened length		Introducer compatibility
			8 ATM nominal pressure	12 ATM rated burst pressure	
12 x 29 mm	85370	85379	25.6 mm	25.3 mm	9 Fr
12 x 41 mm	85371	85380	37.6 mm	37.1 mm	9 Fr
12 x 61 mm	85372	85381	58.2 mm	57.6 mm	9 Fr

The Advanta V12 stent is not available in the U.S. Getinge, **GETINGE**  are trademarks or registered trademarks of Getinge AB, its subsidiaries, or affiliates in the United States or other countries. Copyright 2023 Atrium Medical Corp. All rights reserved. 04/2023 · PN 011529 Rev AB

 **Manufacturer** · Atrium Medical Corporation · 40 Continental Blvd · Merrimack, NH 03054 · USA · +1 603-880-1433

[www.getinge.com](http://www.getinge.com)

**GETINGE** 