

**SAFER AND MORE EFFICIENT
ENDOSCOPE REPROCESSING
WITH THE NEW GETINGE ED-FLOW AER**

Getinge has announced the introduction of the world's first pass-through AER, with two separate, asynchronous chambers. The new model will also enable high-volume endoscopy units to benefit from Getinge's technology for safer and faster reprocessing of flexible endoscopes.

The double-door pass-through unit is designed to be built into a wall, thus assuring a physical barrier soiled and clean areas, thereby reducing the risk of cross-contamination in larger hospitals and endoscopy clinics.

Having two asynchronous chambers and a 22-25-minute process means fast turnaround times and high throughput and up to 6 endoscopes can be run in just above an hour.

The Getinge ED-FLOW AER also offers hands-free operation and RFID technology, which means that it can be loaded, started

Getinge Infection Control
P O Box 69
SE 310 44 Getinge
Sweden
Press Contact:
Susanne Werneman
Phone: +46 35 15 55 00
Fax: +46 35 1814 50
susanne.werneman@getinge.com

and unloaded without the operator having to touch any panels or buttons with soiled hands.

The Getinge ED-FLOW AER complies with EN ISO 15583, parts 1 and 4, and has been tested by the independent test institute Biotech Germande in Marseille.

“With its true barrier system, combined with the flexibility of the two chambers, the Getinge ED-FLOW AER brings a whole new dimension to endoscope reprocessing departments, in terms of both infection control perspective and productivity,” says Jessika Broström, Marketing Manager at Getinge.

The Getinge ED-FLOW AER will be launched at the Medica exhibition in Düsseldorf, 19-22 November 2008. For more information, please contact Marketing Manager, Jessika Broström. Office: + 46 35 15 55 36, Mobile: + 46 703 18 61 62 or E-mail: jessika.brostrom@getinge.com

Getinge Infection Control
P O Box 69
SE 310 44 Getinge
Sweden
Press Contact:
Susanne Werneman
Phone: +46 35 15 55 00
Fax: +46 35 1814 50
susanne.werneman@getinge.com

ABOUT GETINGE INFECTION CONTROL

Getinge is a world leader in complete solutions for effective and efficient cleaning, disinfection and sterilization. Our know-how comprises everything from architectural planning, production and handling equipment, to systems for full traceability of sterile goods. Our commitment covers expert advice, training and long-term technical support.

Encl. Photos.

Getinge Infection Control
P O Box 69
SE 310 44 Getinge
Sweden
Press Contact:
Susanne Werneman
Phone: +46 35 15 55 00
Fax: +46 35 1814 50
susanne.werneman@getinge.com