

New Getinge releases on display atACHEMA

Getinge will be exhibiting three new product atACHEMA, the 29th International Exhibition Congress on Chemical Engineering, Environmental Protection and Biotechnology. This event will be held 11-15 May 2009 in Frankfurt am Main, Germany.

The new products include the Getinge GEV ventilation sterilizer, the STERBOX E-Beam sterilization system, and the new ISOTEST for fast and efficient sterility testing.

The Getinge GEV ventilator sterilizer is specifically designed for sterilization of sterile injectable products, and features a gentle process that also dries and cools products thoroughly.

The unique and innovative E-Beam sterilization system, the STERBOX, is designed for in-line terminal sterilization of small-volume parenterals and medical devices. This new product – and the company's partnership with Compliant Logistics – enables Getinge to provide complete, turnkey solutions for handling parenteral products.

Getinge's brand new isolator, ISOTEST, is being shown publicly for the first time atACHEMA. The system features a completely new integrated HPV biodecontamination system, Steritrace II, for rapid and efficient sterility testing. The new ISOTEST represents a new sterile transfer concept, and can be also adapted for traditional compendial sterility testing or rapid microbiology methods.

Getinge also provides contamination prevention systems for QC/QA (quality control / quality assurance) laboratories, as well as a wide range of laboratory glassware washers and laboratory sterilizers, which will also be on display at the Getinge booth in hall 6.2 (H25-26).

Press contact

Susanne Werneman, Phone +46 35 15 55 00
susanne.werneman@getinge.com

About Getinge Infection Control

Getinge provides complete solutions for effective and efficient cleaning, disinfection and sterilization. The company's know-how comprises everything from architectural planning, production and handling equipment, to systems for full traceability of sterile goods. Getinge's commitment covers expert advice, training and long-term technical support.

www.getinge.com