


# Advanta V12 Large Diameter



Trusted. Reliable. Proven.

## Now Larger.

- Increasing your treatment options with 12mm diameter balloon expandable covered stents.
- Full PTFE encapsulation provides smooth flow lumen.
- Over 15 years of consistent, durable performance.

# Ordering Information

## Advanta V12 balloon expandable covered stent 5 - 10 mm Diameter, .035" guidewire

Stent Diameter/Length	Order Number 80 cm Catheter Length	Order Number 120 cm Catheter Length	Foreshortened Length		Introducer Compatibility
			8 ATM Nominal Pressure	12 ATM Rated Burst Pressure	
5 x 16 mm	85340	85350	15.9 mm	15.6 mm	6 Fr
5 x 22 mm	85341	85351	21.3 mm	21.0 mm	6 Fr
5 x 32 mm	85388	85394	32.3 mm	32.3 mm	7 Fr
5 x 38 mm	85320	85330	37.2 mm	37.7 mm	7 Fr
5 x 59 mm	85321	85331	58.6 mm	60.0 mm	7 Fr
6 x 16 mm	85342	85352	15.7 mm	15.1 mm	6 Fr
6 x 22 mm	85343	85353	20.8 mm	20.2 mm	6 Fr
6 x 32 mm	85389	85395	31.7 mm	31.5 mm	7 Fr
6 x 38 mm	85322	85332	36.6 mm	37.0 mm	7 Fr
6 x 59 mm	85323	85333	57.8 mm	58.7 mm	7 Fr
7 x 16 mm	85344	85354	15.0 mm	14.2 mm	7 Fr
7 x 22 mm	85345	85355	20.1 mm	19.4 mm	7 Fr
7 x 32 mm	85390	85396	31.3 mm	31.2 mm	7 Fr
7 x 38 mm	85324	85334	35.8 mm	35.7 mm	7 Fr
7 x 59 mm	85325	85335	57.1 mm	57.5 mm	7 Fr
8 x 32 mm	85391	85397	30.0 mm	29.6 mm	7 Fr
8 x 38 mm	85326	85336	34.7 mm	34.7 mm	7 Fr
8 x 59 mm	85327	85337	56.0 mm	56.5 mm	7 Fr
9 x 32 mm	85392	85398	28.7 mm	29.2 mm	7 Fr
9 x 38 mm	85328	85338	33.7 mm	32.7 mm	7 Fr
9 x 59 mm	85329	85339	54.6 mm	54.0 mm	7 Fr
10 x 38 mm	85360	85364	30.8 mm	30.9 mm	7 Fr
10 x 59 mm	85361	85365	53.3 mm	52.5 mm	7 Fr

## 12 mm Large Diameter, .035" guidewire

**New Size Available**

Stent Diameter/Length	Order Number 80 cm Catheter Length	Order Number 120 cm Catheter Length	Foreshortened Length		Introducer Compatibility
			8 ATM Nominal Pressure	10 ATM Rated Burst Pressure	
12 x 29 mm	85370	85379	25.6 mm	25.3 mm	9 Fr
12 x 41 mm	85371	85380	37.6 mm	37.1 mm	9 Fr
12 x 61 mm	85372	85381	58.2 mm	57.6 mm	9 Fr

The Advanta V12 balloon expandable covered stent is CE marked and TGA approved for restoring and improving the patency of iliac and renal arteries. Renal approval is for 5-7 mm diameter sizes. Advanta V12 has Canadian Health Ministry license for restoring the patency of iliac lesions. The Advanta V12 stent is not available in the U.S.

The Advanta V12 balloon expandable covered stent is manufactured by Atrium Medical Corporation • Protected by the following international and U.S. patent(s): <http://patents.maquet.com>. • Getinge, **GETINGE** and Atrium are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates in the United States or other countries • Atrium is registered with the U.S. Patent and Trademark Office. • Copyright 2020 Atrium Medical Corp. • All rights not expressly granted are reserved. • Refer to Instructions for Use for current indications, warnings, contraindications, and precautions. • 02/20

Getinge • 40 Continental Blvd • Merrimack, NH 03054 • USA • +1 603-880-1433

PN 011529 Rev AA