

Getinge Assured

Wash Monitor

The Getinge Assured Wash Monitor is a single use synthetic soil test for routine monitoring of automated wash processes.

The test is inserted into a Wash Monitor Holder and placed into the washer disinfectant.



Performance Characteristics

The synthetic test soil (red diamond) simulates human blood and tissue. When inserted into the Wash Monitor Holder the test soil is half obscured, thus simulating the joint of a surgical instrument. The removal of all test soil indicates effective cleaning.

Physical Characteristics

The Wash Monitor consists of a metal strip with synthetic test soil.

- **Dimensions:** 89 mm x 25 mm (3.5" X 1")
- **Substrate Material:** Metal

Environmental Limitations

Storage: 4° - 30°C (40° - 86°F) and 10-60% Relative Humidity. Keep away from disinfectants and chemical fumes.

Shelf Life

24 months. The expiry date is printed on the product packaging.

Standards Classification

Consistent with ISO 15883, routine monitoring of washer disinfectors.

Ordering Information

Getinge Assured	Ordering number
Wash Monitor 100 pcs	503864700
Wash Monitor Holder	503864800