

# Getinge Assured

## Wash Monitor Ultrasonic

The Getinge Assured Ultrasonic Wash Monitor is a single use synthetic soil test for routine monitoring of an ultrasonic wash process.

The test is inserted in a Wash Monitor Holder and then placed in the ultrasonic washer.



### Performance Characteristics

The synthetic test soil (blue diamond) simulates human blood and tissue. The removal of all test soil indicates effective cleaning.

### Physical Characteristics

The Wash Monitor Ultrasonic consists of a metal strip with synthetic test soil.

- **Dimensions:** 89 mm x 25 mm (3.5" X 1")
- **Substrate Material:** Metal

### Environmental Limitations

Storage: 4° - 30°C (40° - 86°F) and 10-60% Relative Humidity. Keep away from disinfectants and chemical fumes.

### Shelf Life

24 months. The expiry date is printed on the product packaging.

### Standards Classification

Consistent with ISO 15883, routine monitoring of washer disinfectors.

### Ordering Information

Getinge Assured	Ordering number
Wash Monitor Ultrasonic 50 pcs	503865200
Wash Monitor Holder	503864800