

# Getinge Assured

## Emulating Indicator 3.5 -A

The Getinge Assured Emulating Indicator 3.5 -A is designed for use in steam sterilizers operating at either 121°C or 134°C.

The Emulating Indicator gives visible indication that sterilizing parameters (time, temperature and steam) were met according to EN ISO 11140-1. The indicator has an adhesive backing for documentation use.



### Performance Characteristics

The Getinge Assured Emulating Indicator 3.5 -A diamond symbol changes color from purple to green when exposed to appropriate steam sterilization conditions.

### Stated Values

Time at 134°C	Time at 121°C
3.5 minutes	9 minutes

\*As determined in a steam sterilization resistometer. Resistometer end point times may be different than standard steam sterilizer endpoint times.

### Physical Characteristics

- **Lamination:** Top
- **Dimensions:** L 38mm x W 38mm
- **Packaging:** 250 indicators in a zip lock moisture barrier bag

### Environmental Limitations

Storage: In provided sealed packaging at room temperature 10°-38°C and 10-60% Relative Humidity. Keep away from light and heat exposure.

### Shelf Life

24 months. The expiry date is printed on the product packaging.

### Standards Classification

EN ISO 11140-1, Type 6.

### Ordering Information

**Getinge Assured**  
Emulating Indicator 3.5 -A

**Order Number**  
504051100