



Cardiosave IABP

Cardiac Support you can count on



Cardiosave

Engineered for Reliability

When patients require hemodynamic support, every moment counts. Whether stabilizing post-AMI, managing cardiogenic shock, heart failure or supporting high-risk PCI, Cardiosave empowers your team with precision, speed, and consistent performance. Built for versatility, it is designed to support patients within the hospital and during transport. With Cardiosave, cardiac support is available when and where it is needed. Because behind every critical moment is a patient.



Cardiosave IABP

Blood Surveillance Algorithm

Monitors for balloon perforation

Automatically Calibrates (zeros)

When using a Getinge Fiber-Optic IAB catheter, Cardiosave automatically calibrates inside the patient on Startup and every 2 hours

IPX1 Rated

Protects electronics from dripping fluid

Color coded keyed connections

to assist setup

Hot-Swappable Batteries

No downtime during AC power loss or transport

Retractable Power Cord

Keeps pathways clear and safe



Focused on Safety and Performance

- Blood surveillance algorithm monitors the IAB catheter to detect the presence of balloon perforation, improving patient safety and reducing risk of damage to the IABP.
- Adaptive automation continuously adjusts and synchronizes timing with dynamic patient physiology.
- Design enhancements to internal components boost overall performance and serviceability.

Designed with the User in Mind

- Hybrid model supports hospital use and converts to Rescue configuration for transport.
- Ergonomic connection panel enables intuitive cable and tubing setup and management in care settings.
- Data integration interfaces with hospital systems to enable EMR integration and real-time monitoring of patient data.



Ordering Information

Product		Part Number
Cardiosave Hybrid - AS/NZS3112 Type I Plug		0998-00-0800-31
Cardiosave Hybrid - Type J Plug		0998-00-0800-32
Cardiosave Hybrid - Type D Plug		0998-00-0800-33
Cardiosave Hybrid - Type K Plug		0998-00-0800-34
Cardiosave Hybrid - Type M Plug		0998-00-0800-35
Cardiosave Hybrid - Type B Plug		0998-00-0800-65 (Japan)

Product		Part Number
Cardiosave Hybrid - Type I Plug		0998-00-0800-45
Cardiosave Hybrid - Type G Plug		0998-00-0800-52
Cardiosave Hybrid - Type B Plug		0998-00-0800-53
Cardiosave Hybrid - Type E/F Plug		0998-00-0800-55
Cardiosave Hybrid - Type N Plug		0998-00-0800-36




**THIS PAGE IS THE INSIDE BACK COVER AND HAS
BEEN LEFT INTENTIONALLY BLANK.**



MCA00004640 Rev A · **GETINGE** ✱, Getinge, and Maquet, are trademarks or registered trademarks of Getinge AB, its subsidiaries, or affiliates in the United States or other countries · Copyright 2026 Datascope Corp. · All rights reserved.
△Refer to Instructions for Use for current indications, warnings, contraindications and precautions · 01/2026

This document is intended to provide information to an international audience outside of the US.

Getinge · 1300 MacArthur Blvd., Mahwah, NJ 07430, USA

 **Manufacturer:** Datascope Corp. 1300 MacArthur Blvd., Mahwah, NJ 07430, USA

www.getinge.com