

Cross-Checks HP

CI 115

The Cross-Checks HP gives visible confirmation that the goods have been exposed to hydrogen peroxide.



Key Features

- Process indicator for monitoring H₂O₂ sterilization cycles
- Indicator can be used as a permanent record

Performance Characteristics

When exposed to H₂O₂ sterilization process, the Cross-Checks HP Indicator will change from purple to a shade of pink.

Physical Characteristics

Dimensions: L 97mm x W 14mm

Color Change

Purple to pink

Storage

The storage conditions are printed on the product packaging. Store in resealable bag provided away from oxidizing agents

Shelf Life

24 months from date of manufacture. The expiry date is printed on the product packaging

SteriTec Part No.

CI 115, Cross-Checks HP

Packaging

24 bags of 250