


Getinge Assured

Green Tape

The Getinge Assured Green Tape can be used in 121°C gravity sterilizers and 134°C pre-vacuum sterilizers.

Stripes turn black after exposure to steam sterilization conditions, thus providing identification of processed items.



Performance Characteristics

The process indicators on the tape will turn from green to black when exposed to steam sterilization.

Physical Characteristics

Green Tape 18mm

- **Dimensions:** L 56m x W 18mm
- **Packaging:** 48 rolls

Green Tape 24mm

- **Dimensions:** L 56m x W 24mm
- **Packaging:** 36 rolls

Environmental Limitations

Storage: At room temperature 10°-38°C and 10-60% Relative Humidity.

Shelf Life

36 months. The expiry date is printed on the product packaging.

Standards Classification

EN ISO 11140-1, Type 1.

Ordering Information

Getinge Assured

Order Number

Green Tape 18mm

6001154902

Green Tape 24mm

6001154903