

Getinge Clean Instrument Lubricant Plus

For both manual and automated lubrication

Use in washer-disinfectors and manual bath applications

Getinge Clean Instrument Lubricant Plus is formulated to efficiently lubricate surgical instrument hinges, pins, joints and box locks for optimal functionality. It also prevents friction corrosion of these moving parts prolonging the life of the instrument.

Instrument Lubricant Plus is silicone and mineral oil free and will not interfere with subsequent steam or gas sterilization processes.



Available in 5L, 25L and 205L containers

This document is intended to provide information to an international audience outside of the US.

For more detailed information, please contact your local Getinge representative at getinge.com/contact

www.getinge.com

GETINGE 