

# Getinge Clean Chamber Foam

— Ready-to-use acid foam spray

Getinge Clean Chamber Foam makes caring for your equipment easier

Getinge Clean Chamber Foam is a ready-to-use acid foam spray that is applied inside process chambers. The foam makes cleaning less time-consuming and is easily wiped off with a wet cloth. This product leaves the chamber clean and free of limescale, mineral deposits and odors. Even rust can be removed.

Getinge Clean Chamber Foam allows for a safer and minimal hands-on cleaning procedure for service technicians and CSSD staff. This product promotes desirable process conditions while minimizing repair costs due to poor water quality and neglected maintenance.

Chamber Foam has been tested on Getinge sterilizers and washer-disinfectors and is compatible with stainless steel surfaces.



---

This document is intended to provide information to an international audience outside of the US.

For more detailed information, please contact your local Getinge representative at [getinge.com/contact](https://www.getinge.com/contact)

[www.getinge.com](https://www.getinge.com)

**GETINGE**