

Maquet Equipment

Conveniently position high-resolution data where you want

During a surgical procedure nothing should disturb the workflow in the operating room. We ensured that the flat screen holder works seamlessly with surgical lights, ceiling service units and operating tables for coordinated performance of staff and equipment. The flat screen holder accommodates larger and heavier flat screens – up to 32 inches, with full HD or 4K visualization, without disrupting your surgical light configuration.



Improving your working environment

- **Safety through clarity**
Bring 4K flat screen imagery to the surgical environment, giving surgeon clear intraoperative visualization of patient tissues.
- **Application-oriented positioning**
Tilt and adjust the flat screens vertically or horizontally to best meet the positioning requirements of each procedure, whether the surgeon is seated or standing.
- **Prevent interference to improve workflow**
The Flat Screen Holder can be ideally positioned during the procedure to not disturb staff or be pushed aside, creating additional space to improve circulating staff activity when screens are not needed.
- **Full coverage of the operating area**
Mounted on Maquet Satellite tube, the screen holder might be placed anywhere within a 4.6 m diameter area centered around the mounting structure.
- **Minimize the risk of hospital-acquired infections**
Reduce the risk of cross-contamination. The hidden cable closed design and dedicated transformer box improve hygiene and minimize the risk of accidental disconnection.



The parking position enhances head clearance in the operating room.



Angulation from + 90° to - 45° to cater surgeon's needs.



The rear box facilitates cleaning on the operating room as it can hold insulation transformers, cables and wireless transmitters.



	Single Basic	Double Basic	Single Premium	Double Premium
Screen compatibility	Up to 32"	Up to 32"	Up to 32"	Up to 32"
Screen holder tilting	N/A	N/A	+90°/-45°	+60°/-10°
Handle	Optional	Optional	Optional	Optional
Rear box	Optional	Optional	Optional	N/A

This information is intended for an international audience outside the US. This information is aimed exclusively at healthcare professionals or other professional audiences and are for informational purposes only, is not exhaustive and therefore should not be relied upon as a replacement of the Instructions for Use, service manual or medical advice. Getinge shall bear no responsibility or liability for any action or omission of any party based upon this material, and reliance is solely at the user's risk. Any therapy, solution or product mentioned might not be available or allowed in your country. Information may not be copied or used, in whole or in part, without written permission by Getinge.

Manufacturer · Maquet S.A.S · Parc de Limère · Avenue de la Pomme de Pin · CS 10008 Ardon · 45074 Orléans, cedex 2 · France
 © 2026 Getinge | Getinge and **GETINGE** * are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates.
 DMS-0019277 v1 | All rights reserved.

