

Getinge Clean

Renuzyme Ultra

Renuzyme Ultra is a high performance enzyme detergent specifically formulated for use at lower temperatures (65-140°F). It combines next generation lipase and protease enzymes to rapidly digest a broad spectrum of protein and lipid based soils. The novel surfactant package provides superior wetting to clean the most intricate of device assemblies. Its highly concentrated formula means economical cost per use in any application.

A Getinge Green Product



Application

Renuzyme Ultra is designed for use in manual sink applications, ultrasonic cleaners and washer-disinfectors. Its highly concentrated formula means economical use with little foaming, ensuring easy visibility in sinks and non-cavitation of pumps in washer-disinfectors.

Key Features

- Formulated for use at lower temperatures
- Concentrated Formula - economical cost per use
- Superior digestion of proteinaceous and lipidic soils
- Low foaming and non-clouding
- Getinge Green Product

Usage Information

Please consult Instructions for Use (IFU) for dilution rates and specific usage information. Equipment manufacturer IFU should be consulted for cleaning instructions.

- Do not add bleach, detergents, or any other chemical to this product. Added chemical may destroy the effectiveness of this cleaner.

Quality Statement: Confidence in Getinge is the most important quality criteria. This is the hallmark of all our external and internal commitments, activities and products. Products and services supplied by Getinge conform to the agreed terms and expectations. The achievement of these quality goals is the basis for continued competitive and successful enterprise.

Physical and Chemical Properties

- No dyes added - clear, slight yellow to amber liquid
- No fragrance added
- pH of concentrate: 7–8.5
- Viscosity similar to water

NOTICE: For institutional and industrial use only. Not for home use. Refer to Safety Data Sheet (SDS) for more information regarding safety, handling, storage and waste disposal methods.

Ordering Information

Getinge Clean	Order Number
4 x 1 Gallon	61301606716
2.5 Gallons	61301606791
5 Gallons	61301606717

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2018 Getinge AB or its subsidiaries or affiliates • All rights reserved. • MCV00044132 REV B

Manufacturer: Intercon Chemical Co. 1100 Central Industrial Dr., St. Louis, MO 63110
Sales Office: Getinge USA Sales, LLC, 1777 East Henrietta Road, Rochester, NY 14623

www.getinge.com

GETINGE 