


Accessories

For Getinge Thoracic Chest Drainage Systems



19910 — Red Cap

- For capping patient line before disposal
- Twenty-five (25) each per case — sterile



19915 — Pediatric Connector 1/4 inch x 1/4 inch

- Larger connector included with infant/pediatric drain
- For use with an adult drain and smaller catheters
 - 16 French or larger thoracic catheters
- Twenty-five (25) each per case — sterile



19918 — Step Down Connector 3/16 inch x 3/8 inch

- For use with adult chest drain
- Fits smaller catheters — 8, 10, 12 French
- Twenty-five (25) each per case — sterile

None of our drain accessories are made with natural rubber latex.

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2025 Getinge AB or its subsidiaries or affiliates • All rights reserved • ⚠ CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician. Refer to Instructions for Use for current indications, warnings, contraindications, and precautions.

Sales Office • Getinge USA Sales LLC • 45 Barbour Pond Dr • Wayne, NJ 07470

www.getinge.com