

Accessories

For Getinge Thoracic Chest Drainage Systems

From admission to discharge, Getinge's portfolio of Atrium brand thoracic drainage products provides all of your therapeutic device needs. Providing everything for initiation, collection, and transition for the changing needs of the patient during therapeutic treatment in the ER, OR, ICU, at bedside, and in transport.

Having the flexibility to change the drainage therapy quickly is key to improving patient outcomes. Atrium thoracic drainage accessories provide a variety of options: catheter-to-chest drain, suction-to-chest drain configurations, and step-down and capping of patient line before disposal.




19910 — Red Cap

- For capping patient line before disposal
- For use if only one patient line is needed on a dual drain
- Twenty-five (25) each per case — sterile



19913 — Adult/Pediatric Connector 1/4 inch x 3/8 inch

- For use with an adult drain and smaller catheters — 16, 20 French
- Twenty-five (25) each per case — sterile



19914 — Standard Adult Connector 3/8 inch x 3/8 inch

- Standard 5-in-1 connector
- For use on larger catheters — 24 French and larger
- Twenty-five (25) each per case — sterile



19915 — Pediatric Connector 1/4 inch x 1/4 inch

- Larger connector included with infant/pediatric drain
- For use with an adult drain and smaller catheters
 - 16 French or larger thoracic catheters, 9 and 11 mm mediastinal catheters
- Twenty-five (25) each per case — sterile



19916 — Sims Connector 3/16 inch x 1/4 inch

- Same connector included with infant/pediatric drain
- For use with smaller catheter — 8, 10, 12 French
- Twenty-five (25) each per case — sterile

Accessories: For Getinge Thoracic Chest Drainage Systems



19918 — Step Down Connector 3/16 inch x 3/8 inch

- For use with adult chest drain
 - Fits smaller catheters — 8, 10, 12 French
 - Twenty-five (25) each per case — sterile
-



19907 — Patient Tube Pack

- 3/8 inch ID tubing
 - Safety stop if patient tube sterility is accidentally compromised
 - Provides cost savings to patient and hospital
 - Not made with natural rubber latex
 - Ten (10) each per case — sterile
-



19927 — Luer-Lock Connector with Tube

- Allows connection of a standard chest drain or mobile drain to a luer-lock catheter
 - Twenty-five (25) each per case — sterile
-



19928 — Suction Y-Adapter

- For use when connecting 2 chest drains to 1 suction source
 - Ten (10) each per case — sterile
-



19930 – Express™ Suction Bulb

- For use with the Express™ Mini 500 (No. 16400)
 - Allows portable, manual vacuum application
 - Ten (10) each per case — sterile
-

All of these accessories are not made with natural rubber latex

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries. · Copyright 2020 Getinge AB or its subsidiaries or affiliates · All rights reserved · ⚠ CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician. Refer to Instructions for Use for current indications, warnings, contraindications, and precautions.

Getinge · 1 Geoffrey Way · Wayne, NJ 07470 · Phone: (888) 9GETUSA
Sales Office, Canada · 90 Matheson Blvd West, Suite 300 · Mississauga, Ontario L5R 3R3 · Canada

www.getinge.com

GETINGE 

MCV00052245 REV B
EXP11/22