

GET2060

Blanket Warming Cabinet

A single chamber blanket warming cabinet with two tempered glass doors. The blanket warmer is warmed using a low-heat density electrothermal cable array. The electrothermal cable is positioned in the floor and two sides of the blanket warmer, providing even heating of the interior chamber. The blanket chamber has an adjustable temperature range (37 °C to 93 °C/98 °F to 200 °F).



Application

The Getinge GET2060 is designed to elevate blanket temperatures to a level which increases patient comfort.

Construction

- Stainless Steel exterior and interior
- Two doors with handle and hinges
- Two epoxy-coated blanket support assemblies
- Right-hinged is standard
- Casters included

Dimensions

Model	Exterior (H x W x D)	Interior (H x W x D)	Net Weight
GET2060	71"(1803mm) x 32"(813mm) x 29.25"(743mm)	56.75"(1441mm) x 26.5"(673mm) x 23.75"(603mm)	347 lbs (157 kg)

Capacity

Model	Bottles	Bags	Blankets
GET2060	N/A	N/A	60-70

Installation

GET2060 is a stand-alone warming cabinet.

Number of doors

- Double door with windows

Installation alternatives

- Foot Kit
- Left-hinged
- Solid Door
- Pass-Through Design

Other installation options

- Manual lock
- Seismic Mounting brackets (specific installation requirements apply)
- Built-in Trim Kit

Clearance Requirements

- 4" (102mm) from rear
- 4" (102mm) from top
- 2" (51mm) from sides

Safety features

- The controller stores all settings and resumes operation using these settings in the event of a power failure.
- A warming shut-off system, separate from the controller, prevents overheating.

Other features

- Integrated controller lock-out
- Timer
- Operate in Celsius or Fahrenheit
- LED display

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2019

Getinge AB or its subsidiaries or affiliates • All rights reserved • MCV00076068 REV B • EXP 02/21

Sales Office: Getinge USA Sales LLC., 1777 East Henrietta Road, Rochester, NY 14623

Manufactured By: Enthemics Medical Systems, W164 N9221 Water Street, Menomonee Falls, WI 53051

www.getinge.com

GETINGE 