

Getinge Clean

PURE® Substrate Compatibility Matrix

	Aluminum	Anodized Aluminum	Titanium	Nickel	Brass	Copper	Matte Chromium Plated
PURE	C	C	A	A	D	D	A

	Glass	Ceramic	Stainless Steel	Buna-n	Viton	PTFE	Polypropylene
PURE	A	A	A	A	A	A	A

	HDPE	Nylon	Polycarbonate	EPDM	Ultem	Radel	Delrin
PURE	A	A	A	A	A	A	B

- A — Excellent
- B — Good: Minor Effect, slight corrosion, or discoloration.
- C — Fair: Moderate Effect, not recommended for continuous use.
Softening or loss of strength, and swelling may occur.
- D — Severe Effect: Not recommended for any use.

Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. • PURE Hard Surface is a registered trademark of PURE Bioscience, Inc. • Copyright 2020 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Sales Office, US • 45 Barbour Pond Drive • Wayne, NJ 07470
 Sales Office, Canada • 90 Matheson Blvd West, Suite 300 • Mississauga, Ontario L5R 3R3 • Canada

www.getinge.com

MCV00105428 REVA
EXP05/22