

Maquet Moduevo Sky Distributor

Optimizing space

After the Slim, Mezzo and Twin distributors, Getinge's range of ceiling supply units, Maquet Moduevo is adding its portfolio with a horizontal solution: Sky.

In Operating Rooms, Sky can be used for anesthesia stations and, thanks to its optimized head clearance, the distributor can be positioned in a safe way.

Available in different versions, Sky can also be equipped with adjustable shelves in order to support medical devices as well as the complete Maquet Moduevo range of accessories.



Maquet Moduevo Sky's benefits include:

Workflow improvements

- Enlarging workplaces: thanks to its optimized head clearance, Maquet Moduevo Sky reduces footprint and allows medical staff to work in a larger environment.
- Usability and ergonomics: outlets and sockets directly available underneath the module.

Covering all needs

- Maquet Moduevo Sky is available in three versions: Classic, Advanced and Advanced with racks. All of them can offer up to three choices of panel lengths to suit your needs (18.5", 26.4" and 34.3"/ 470, 670 and 870 mm).

Wide choice of options

- An optional underneath mounted ambient light creates a soothing atmosphere, facilitating the work and recording of the medical staff.
- Compatible with all Maquet Moduevo accessories thanks to the EPoS interface: handles, shelves, infusions equipment holders, cable management solutions.
- Shelf loading up to 110.23 lbs./50 kg adjustable height.
- Mechanical bars and EPoS interfaces on racks (both sides).
- Large shelf with control function.
- Smart accessory cable holder for cable management.



Three panels for outlets and sockets (including underneath).



Full range of accessories available.



Smart cable management.



Sky Classic version



Sky Advanced version



Sky Advanced version with racks

Models	Sky Classic version	Sky Advanced version	Sky Advanced version with racks
Length	18.5", 26.4" and 34.3" (470, 670 and 870 mm)		18.5", 26.4" and 34.3" (470, 670 and 870 mm) with 15.7" (400 mm) racks 26.4" and 34.3" (670 and 870 mm) with 31.5" (800 mm) racks (two shelves max.)
Panels for sockets & outlets		One underside and two big side panels (rear and front)	
Accessories positioning	2 lateral mechanical bars (no EPoS) 132.27 lbs./60 kg net load each	2 EPoS: 132.27 lbs./60 kg net load each	2 EPoS: 264.55 lbs./120 kg net load in global
Beams compatibility	Fixed and rotary columns Mechanical beams (without electro-pneumatic or electro-magnetic brakes)	Every Maquet Moduevo ceiling pendant (mechanical beams and motorized beam: Energy) with or without brakes	

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2020 Getinge AB or its subsidiaries or affiliates • All rights reserved - ⚠ CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician. Refer to Instructions for Use for current indications, warnings, contraindications, and precautions.