

# Heli-Check FORM

## LT 112

The Heli-Check FORM verifies air removal, formaldehyde exposure and deep vacuum achievement necessary for the sterilization of hollow lumen instruments.

The Heli-Check FORM is suitable for daily sterilizer monitoring. The test is designed to be used in each sterilization cycle as an independent control device.



### Key Features

- Reusable device
- Adhesive backing for permanent record
- Lead and toxic heavy metal free

### Performance Characteristics

The Heli-Check FORM is a reusable device accompanied with single use chemical indicator strips. The indicator strips will change color from red to green when exposed to appropriate sterilization conditions. The adhesive indicator strips can be used for record keeping.

### Physical Characteristics

Dimensions strip: L 77mm x W 6mm

### Color Change

Red to green for complete change

### Storage

The storage conditions are printed on the product packaging. Store in the re-sealable bag provided away from oxidizing agents

### Shelf Life

24 months from date of manufacture. The expiry date is printed on the product packaging

### SteriTec Part No.

LT 112, Heli-Check FORM

### Packaging

6 bags of 250 tests & 1 device