

# Wash-Checks E Holder

## WC 107

The Wash-Checks E Holder is an auxiliary device for routine monitoring of automated wash processes for hollow instruments in washer-disinfectors and ultrasonic washers.

Wash-Checks E Holder should only be used with Wash-Checks H indicators (in washer-disinfectors) or Wash-Checks UH indicators (in ultrasonic washers). The Wash-Checks H or UH is inserted into the Wash-Checks E Holder, simulating contaminated flexible hollow lumen instruments.



### Key Features

- Reusable holder
- Electro polished stainless steel AISI 304
- Silicone conforms to European specification ST-EC-60-001
- Use with WC 105 and WC 110

### Performance Characteristics

The Wash-Checks H or Wash-Checks UH is inserted into the Wash-Checks E Holder. The holder is then attached to a luer lock connector located on the washer rack.

### Physical Characteristics

- **Dimensions:** 1067 mm X 11 mm
- **Weight:** 73 g

### Shelf Life

Not applicable

### SteriTec Part No.

WC 107, Wash-Checks E Holder

### Packaging

Each