

Customer Notification Letter

Please forward this information to all relevant users and biomedical staff at your facility.

Longer User interface cable for Servo-u and Servo-n ventilator systems

Dear Customer,

In view of the situation with the Coronavirus disease COVID-19, and your institutes recent inquiry, we have designed a longer User interface cable between the patient unit and the User interface of the Servo-u and Servo-n ventilator systems.

The cable delivered as standard with the Servo-u and Servo-n is 2.9 meter / 9.5 feet long.

The new optional cable is 6 meter / 20 feet long. This will allow a wider separation of the User interface from the Patient unit to further enhance control of the ventilator by clinicians outside of the highest risk area during mechanical ventilation of COVID-19 patients.

| Part number | Description |
|-------------|-----------------------------|
| 68 90 333 | Control Cable 6 m Servo-u/n |

In addition, we would like to remind you that the Servo-u/ Servo-n User interface holder is required when mounting of the User interface on a rail 10 x 25-35 mm or a table with 10-45 mm thickness, and must be ordered separately. Its Part number is 68 80 099.

Please be aware that we introduce the alternative cable upon your request for Emergency use purpose during the COVID-19 pandemic and that any purchase will be made under this condition. To order it, please contact your local Getinge representative

Yours Sincerely,

Jennie Haag
Director of Global Product Management
Maquet Critical Care AB

Miray Kärnekull
Medical Director
Maquet Critical Care AB