

Maquet Lyra

A stable, flexible mobile OR table

Sometimes, a reliable OR table is all your hospital needs. The Maquet Lyra Mobile OR Table offers the stability and functionality needed across surgical disciplines, all for a comfortable price.



The advantages of Maquet Lyra include:

Excellent balance of price and performance

Why pay more for features that you don't need? This cost-effective mobile OR table is flexible and easy to use, comfortably positioning patients for most types of surgical procedures.

Exceptional stability and safety

The robust base ensures stability and safety for patients weighing up to 360 kg (794 lbs). The Maquet Lyra Table is capable of a broad range of angles and tilts for extraordinary cross-disciplinary flexibility.

Versatility across surgical disciplines

By combining universal usage features with a broad range of accessories, the Maquet Lyra Mobile OR Table is a practical table for all-around use. Maquet Lyra features the same interfaces as past Getinge Mobile OR tables and the compatibility with existing Getinge accessories ensures seamless integration into your surgical suite.

Reliability and user-friendly operation:

The special design of the castors allows the table to be used immediately, without being lowered onto a baseplate. The free-wheeling castors give the flexibility to move in both straight lines or laterally. Fewer electrical components simplify maintenance for cost-effective long-term use.

Longitudinal shift improves access for better outcomes

The manual longitudinal shift improves intraoperative diagnostics and treatment. It delivers both increased access for the surgeon, and simultaneously improves accessibility for the C-arm. 360° imaging is possible when paired with carbon fiber accessories.

Construction features

- Stable base construction with four double swivel castors for easy movement and maneuvering (base can be locked via foot pedal)
- Base cover made of SMC (sheet moulding composite)
- The override panel is connected to a separate circuit in order to ensure that the table can be safely controlled even in the event of a major software blackout
- Identical interfaces on normal and on reverse side

Electrical specifications

- Specially designed rechargeable batteries, with a single charge lasting about one week in the operating room
- Electronic charge monitoring with visual indicators
- Batteries recharged from the mains supply
- 100–240 V AC (adjustable), 50–60 Hz, over mains cable
- Safety class II, type B; the enclosure leakage current meets the requirements of the patient leakage current for type B conditions as per EN 60601-1

Characteristics Maquet Lyra 1118.10XX

Net weight (without accessories)	190 kg / 419 lbs
Length of OR table without accessories	745 mm / 29"
Maximum overall load (patient plus accessories)	360 kg / 794 lbs with restrictions 180 kg / 397 lbs without restrictions
Trendelenburg / Reverse Trendelenburg	25° / 25°
Lateral tilt	15°
Lower back plate	+75° / -40°
Leg plate (manual)	-90° (down)
Manual longitudinal shift	310 mm / 12"



This document is intended to provide information to an international audience outside of the US.
Maquet (Suzhou) Co., Ltd. · No.158 Fangzhou Road · Suzhou Industrial Park · 215024 Suzhou · China · +86 (0) 512 62839880