

Ordering Information

Advanta V12 balloon expandable covered stent

5-10 mm diameter, .035" guidewire

Stent Diameter (mm)	Stent Length (mm)	Order Number 80 cm catheter length	Order Number 120 cm catheter length	Post Nominal Pressure 8ATM (mm)		Post Rated Burst Pressure 12ATM (mm)		Maximum Recommended Oversizing diameter (mm)	Introducer Compatibility
				Outer Diameter	Length	Outer Diameter	Length		
5	16	85340	85350	4.9	15.9	5.2	15.6	7	6 FR
5	22	85341	85351	4.9	21.3	5.2	21	7	6 FR
5	32	85388	85394	5	32.3	5.3	32.3	9	7 FR
5	38	85320	85330	5.1	37.2	5.3	37.7	9	7 FR
5	59	85321	85331	5	58.6	5.3	60	9	7 FR
6	16	85342	85352	5.7	15.7	6.2	15.1	7	6 FR
6	22	85343	85353	5.8	20.8	6.2	20.2	7	6 FR
6	32	85389	85395	5.9	31.7	6.3	31.5	9	7 FR
6	38	85322	85332	6	36.6	6.3	37	10	7 FR
6	59	85323	85333	6	57.8	6.3	58.7	10	7 FR
7	16	85344	85354	6.9	15	7.3	14.2	7	7 FR
7	22	85345	85355	6.9	20.1	7.3	19.4	7	7 FR
7	32	85390	85396	6.9	31.3	7.3	31.2	9	7 FR
7	38	85324	85334	6.9	35.8	7.3	35.7	10	7 FR
7	59	85325	85335	7	57.1	7.3	57.5	10	7 FR
8	32	85391	85397	7.9	30	8.4	29.6	9	7 FR
8	38	85326	85336	8.1	34.7	8.5	34.7	10	7 FR
8	59	85327	85337	8	56	8.4	56.5	10	7 FR
9	32	85392	85398	8.9	28.7	9.3	29.2	9	7 FR
9	38	85328	85338	8.9	33.7	9.3	32.7	10	7 FR
9	59	85329	85339	8.9	54.6	9.3	54	10	7 FR
10	38	85360	85364	10	30.8	10.4	30.9	10	7 FR*
10	59	85361	85365	9.9	53.3	10.3	52.5	10	7 FR*


12 mm diameter, .035" guidewire

Stent Diameter (mm)	Stent Length (mm)	Order Number 80 cm catheter length	Order Number 120 cm catheter length	Post Nominal Pressure 8ATM (mm)		Post Rated Burst Pressure 10ATM (mm)		Maximum Recommended Oversizing diameter (mm)	Introducer Compatibility
				Outer Diameter	Length	Outer Diameter	Length		
12	29	85370	85379	11.1	25.6	11.4	25.3	16	9 FR
12	41	85371	85380	11.2	37.6	11.6	37.1	16	9 FR
12	61	85372	85381	11.2	58.2	11.6	57.6	16	9 FR

*If using a 10mm device with a Cook introducer sheath, it is recommended to use an 8 French sheath to minimize potential for resistance.

The Advanta V12 stent is not available in the U.S.

The Advanta V12 balloon expandable covered stent is manufactured by Atrium Medical Corporation

• Getinge and **GETINGE**  are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates in the United States or other countries • Copyright 2024 Atrium Medical Corp. • All rights not expressly granted are reserved. • Refer to Instructions for Use for current indications, warnings, contraindications, and precautions. • PN 011529 Rev AC 09/2024

 **Manufacturer:** Atrium Medical Corporation • 40 Continental Blvd • Merrimack, NH 03054 • USA • +1 603-880-1433