

# Ordering Information

## Advanta V12 covered stent system

5-10 mm diameter, .035" guidewire

Stent Diameter (mm)	Stent Length (mm)	Order Number 80 cm catheter length	Order Number 120 cm catheter length	Post Nominal Pressure 8ATM (mm)		Post Rated Burst Pressure 12ATM (mm)		Maximum Recommended Post-Dilation Balloon Diameter (mm) <sup>1</sup>	Introducer Compatibility
				Outer Diameter	Length	Outer Diameter	Length		
5	16	GBX05161C	GBX05162C	4.9	15.9	5.2	15.6	7	6 FR
5	22	GBX05221C	GBX05222C	4.9	21.3	5.2	21	7	6 FR
5	32	GBX05321C	GBX05322C	5	32.3	5.3	32.3	9	7 FR
5	38	GBX05381C	GBX05382C	5.1	37.2	5.3	37.7	9	7 FR
5	59	GBX05591C	GBX05592C	5	58.6	5.3	60	9	7 FR
6	16	GBX06161C	GBX06162C	5.7	15.7	6.2	15.1	7	6 FR
6	22	GBX06221C	GBX06222C	5.8	20.8	6.2	20.2	7	6 FR
6	32	GBX06321C	GBX06322C	5.9	31.7	6.3	31.5	9	7 FR
6	38	GBX06381C	GBX06382C	6	36.6	6.3	37	10	7 FR
6	59	GBX06591C	GBX06592C	6	57.8	6.3	58.7	10	7 FR
7	16	GBX07161C	GBX07162C	6.9	15	7.3	14.2	7	7 FR
7	22	GBX07221C	GBX07222C	6.9	20.1	7.3	19.4	7	7 FR
7	32	GBX07321C	GBX07322C	6.9	31.3	7.3	31.2	9	7 FR
7	38	GBX07381C	GBX07382C	6.9	35.8	7.3	35.7	10	7 FR
7	59	GBX07591C	GBX07592C	7	57.1	7.3	57.5	10	7 FR
8	32	GBX08321C	GBX08322C	7.9	30	8.4	29.6	9	7 FR
8	38	GBX08381C	GBX08382C	8.1	34.7	8.5	34.7	10	7 FR
8	59	GBX08591C	GBX08592C	8	56	8.4	56.5	10	7 FR
9	32	GBX09321C	GBX09322C	8.9	28.7	9.3	29.2	9	7 FR
9	38	GBX09381C	GBX09382C	8.9	33.7	9.3	32.7	10	7 FR
9	59	GBX09591C	GBX09592C	8.9	54.6	9.3	54	10	7 FR
10	38	GBX10381C	GBX10382C	10	30.8	10.4	30.9	10	7 FR <sup>2,3</sup>
10	59	GBX10591C	GBX10592C	9.9	53.3	10.3	52.5	10	7 FR <sup>2,3</sup>

12 mm diameter, .035" guidewire

Stent Diameter (mm)	Stent Length (mm)	Order Number 80 cm catheter length	Order Number 120 cm catheter length	Post Nominal Pressure 8ATM (mm)		Post Rated Burst Pressure 10ATM (mm)		Maximum Recommended Post-Dilation Balloon Diameter (mm) <sup>1</sup>	Introducer Compatibility
				Outer Diameter	Length	Outer Diameter	Length		
12	29	85370	85379	11.1	25.6	11.4	25.3	16 <sup>4</sup>	9 FR
12	41	85371	85380	11.2	37.6	11.6	37.1	16 <sup>4</sup>	9 FR
12	61	85372	85381	11.2	58.2	11.6	57.6	16 <sup>4</sup>	9 FR

<sup>1</sup>Note: For maximum stent diameters and associated foreshortened lengths, see the Advanta V12 covered stent system Instructions for Use.

<sup>2</sup>If using a 10mm device with a Cook introducer sheath, it is recommended to use an 8 French sheath to minimize potential for resistance.

<sup>3</sup>If using a 10 mm device with a Cook introducer sheath through a femoral access for a BEVAR procedure, it is recommended to use a 9 French sheath to minimize potential for resistance.

<sup>4</sup>Stents post-dilated using a 16mm balloon have not been tested in kissing stent configuration.

• Getinge and **GETINGE**  are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates in the United States or other countries • Copyright 2024 Atrium Medical Corporation. • All rights not expressly granted are reserved. • Refer to Instructions for Use for current indications, warnings, contraindications, and precautions. • Advanta V12 covered stent system is not available in the U.S.

PN 011529 Rev AD 12/2025

 **Manufacturer:** Atrium Medical Corporation • 40 Continental Blvd • Merrimack, NH 03054 • USA • +1 603-880-1433

[www.getinge.com](http://www.getinge.com)

**GETINGE** 