




**Getinge  
Vacuum  
Stabilizer**

OM-10000  
OM-9100S  
OM-9000S



**Max. setting  
400 mmHg**

**GETINGE**

Use this setting for these  
vacuum stabilizers:

|             |          |
|-------------|----------|
| Acrobat-i   | OM-10000 |
| Acrobat V   | OM-9100S |
| Acrobat SUV | OM-9000S |

How can we help?

Contact us: 888 627 8383

Sales Office: Getinge  
45 Barbour Pond Drive  
Wayne, NJ 07470 USA  
[www.getinge.com](http://www.getinge.com)

⚠ CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician. Refer to Instructions for Use for current indications, warnings, contraindications, and precautions.

MCV00007129 REV B