




**Getinge
Vacuum
Positioner**

XP-5000
XP-3000



**Max. setting
250 mmHg**

GETINGE 

Use this setting for these
vacuum positioners:

Acrobat-i	XP-5000
Xpose 3	XP-3000

How can we help?
Contact us: 888 627 8383

Sales Office: Getinge
45 Barbour Pond Drive
Wayne, NJ 07470 USA
www.getinge.com

⚠ CAUTION: Federal (US) law restricts this device
to sale by or on the order of a physician. Refer to
Instructions for Use for current indications,
warnings, contraindications, and precautions.

MCV00007495 REVB