

Getinge Clean Chamber Foam

– Ready-to-use acid foam spray

Getinge Clean Chamber Foam makes caring for your autoclave easier.

The Getinge Clean detergent product portfolio has expanded with a new cleaning product: Getinge Clean Chamber Foam, a ready-to-use acid foam spray that is applied inside process chambers. The foam makes cleaning less time-consuming and is easily wiped off with a wet cloth. The product leaves the chamber clean and free of limestone, mineral deposits and odors. Even rust can be removed.

Getinge Clean Chamber Foam allows for a safer and minimum hands-on cleaning procedure for service technicians and CSSD staff. This product will promote desirable process conditions whilst minimizing repair costs due to poor water quality and neglected maintenance.

The product has been tested on Getinge sterilizers and washer-disinfectors and is compatible with all stainless steel surfaces.

This document is intended to provide information to an international audience outside of the US.

www.getinge.com



2916-getinge-clean-chamber-foam-leaflet-a5-170413-en

GETINGE 