


# Maquet Exali 10 Examination Light

## High-performance LED lighting engineered for multiple medical specialties

GETINGE

# Maquet Exali 10 Examination Light



## Performant illumination

- 50,000 lux at 0.5 m / 1.64 ft for precise examinations
- Uniform 19 cm / 7.48 in light field for consistent visibility
- CRI 96 for excellent color accuracy

## Easy & flexible to use

- One-handed positioning
- Gooseneck or spring-balanced arm options
- 300° rotation and 180° tilt
- Up to 1.3 m / 4.27 ft reach with extension arm
- Intuitive touch control

## Engineered to drive down lifetime costs

- 50,000-hour LED lifetime
- No maintenance or bulb replacement
- Energy efficient (~18 W)
- Durable ABS / polymer housing

## Optimized for many specialties

- 3 adjustable color temperatures: 3,300 K / 3,800 K / 4,400K
- 5 adjustable brightness levels

## Designed to prevent cross contamination

- Sealed light head
- Smooth, disinfectant-compatible surfaces

## Flexible installation

- Rail clamp
- Wall mount
- Mobile stand

This information is intended for an international audience outside the US. This information is aimed exclusively at healthcare professionals or other professional audiences and are for informational purposes only, is not exhaustive and therefore should not be relied upon as a replacement of the Instructions for Use, service manual or medical advice. Getinge shall bear no responsibility or liability for any action or omission of any party based upon this material, and reliance is solely at the user's risk. Any therapy, solution or product mentioned might not be available or allowed in your country. Information may not be copied or used, in whole or in part, without written permission by Getinge.

Manufacturer · Derungs Licht AG. Hofmattstrasse 12. CH-9200 Gossau SG. Switzerland

Maquet Exali is a registered trademark.

© 2026 Getinge | Getinge and GETINGE are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates.

DMS-0008785 v4 | All rights reserved.

GETINGE