

Maquet Moduevo Bridge Ceiling Supplys Units

Think Horizontal. Put Care First.

Maquet Moduevo Bridge redefines ceiling supply by bringing critical equipment exactly where caregivers need it — close, visible, and accessible. Its horizontal design maximizes space, preserves full patient access, and enables confident care delivery at every acuity level.

Optimized for:

- Intensive Care Units
- High Dependency Units
- Recovery Areas
- Emergency Department

GETINGE

Maquet Moduevo Bridge

Optimize Space. Elevate Care.

Designed to overcome facility constraints

- Engineered for low-ceiling environments: adjustable beam height from 1800 to 2350 mm (70.9”–90.5”) ensures full functionality without compromising room constraints.
- Built to scale with large units: place multiple bridges side by side and combine beam lengths from 2200 to 3500 mm (86.6”–137.8”)* to maximize capacity without sacrificing valuable space.
- Optimized for windowless rooms: the optional Somnus light* delivers advanced daylight simulation to support patient well-being and staff comfort.
- Designed for evolving care demands: the Ergonomic Positioning System (EPoS) enables fast, tool-free addition or removal of equipment for positioning, holding, and storage — whenever workflows change.

Engineered for Smarter Resource Utilization

- Serve two beds with one unit, maximizing flexibility and optimizing investment during peak demand.

Fostering a healing environment

- Unrestricted access where it matters most: the horizontal design ensures immediate, unobstructed access to the patient’s head, supporting fast, efficient, and confident care delivery.
- Stronger clinical visibility on demand: the optional OR light adapter* enables seamless integration of additional examination lighting using Maquet Ezea or Maquet Volista Access.
- A calmer environment for patients and staff: integrated ambient lighting with up to six color options helps create a soothing atmosphere for patients, families, and caregivers alike.
- Quiet by design: low-noise electromagnetic brakes support a restful care setting, reducing disturbances for patients while maintaining precise positioning.

**Only for Rack and Pole, Value and Premium versions*



The Sommus light has been developed to the Lighting Research Center Recommandations. Bluish tones simulate morning light, while reddish ones represent evening light (based on the Circadian Stimulus).



The lateral slide can be managed by one person. Distributors can be easily moved closer of farther away with minimal effort.



Rack & pole version
A smart and efficient bedside entry solution combining essential gas, power, and lighting in a clean, intuitive design for everyday use.



Q version
A streamlined solution that enhances flexibility and bedside organization, supporting ICUs that require greater capability with fundamental control.



Value version
A versatile ICU solution engineered to improve caregiver workflow, offering enhanced ergonomics, better patient-head access, and integrated examination lighting at the bed side.



Premium version
A premium, fully featured platform designed to deliver maximum ergonomics, advanced Moduevo capabilities, and superior access and flexibility for high-acuity ICUs.

Maquet Moduevo Bridge

Technical data

Technical data	Rack & pole solution	Q version	Value solution	Premium solution
Recommended height clearance	1800-2000 mm (70.8"-78.7")	2100 mm (82.6")	2100 mm (82.6")	2350 mm (92.5")
Beam lenghts	2200, 2500, 2800, 3100, 3500 mm (86.6", 98.4", 110.2", 137.8")	2200, 2500, 2800, 3100 mm (86.6", 98.4", 110.2")	2200, 2500, 2800, 3100, 3500 mm (86.6", 98.4", 110.2", 137.8")	2200, 2500, 2800, 3100, 3500 mm (86.6", 98.4", 110.2", 137.8")
Load capacity	Pole: 60 kg / 132 lbs (30 kg / 66 lbs on each side) Rack: 130 kg / 286 lbs	150 kg / 331 lbs for each shuttle	150 kg / 331 lbs for each shuttle	With 400 mm swivel arm : 120 kg / 264 lbs With 600 mm swivel arm : 80 kg / 176 lbs
Braking systems	Mechanical	Mechanical	Mechanical & Electromagnetic	Mechanical & Electromagnetic
Lighting systems	Top night light Reading or task light Somnus light and colored ambient light	Top night light Reading or task light	Top night light Reading or task light Somnus light and colored ambient light	Top night light Reading or task light Somnus light and colored ambient light
Bridge adaptor for OR Lights Maquet Ezea or Maquet Volista Access II	yes	n/a	yes	yes

This information is intended for an international audience outside the US. This information is aimed exclusively at healthcare professionals or other professional audiences and are for informational purposes only, is not exhaustive and therefore should not be relied upon as a replacement of the Instructions for Use, service manual or medical advice. Getinge shall bear no responsibility or liability for any action or omission of any party based upon this material, and reliance is solely at the user's risk. Any therapy, solution or product mentioned might not be available or allowed in your country. Information may not be copied or used, in whole or in part, without written permission by Getinge.

Manufacturer · Maquet (Suzhou) Co., Ltd. · No.158 Fangzhou Road, Suzhou, P.R. China 215024 · +86 (0) 512 62839880

© 2026 Getinge | Getinge and GETINGE are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates.
DMS-00015775 v1 | All rights reserved.

GETINGE